



Highly reproducible Illumina libraries

Volta's DNA Library Prep for Illumina™ App offers an automated, streamlined process for efficiently converting DNA of variable quantities into Illumina™ compatible sequencing-ready libraries with minimal bias and artifacts. The App supports whole genome sequencing, whole exome sequencing, and targeted sequencing.

Benefits

- Delivers consistent fragment lengths
- Low duplication rates of <10%
- Robust workflow with <1 in 300 failure rate

App specifications

Input: **50 – 600 ng DNA**

Output: **Sequencing-ready dsDNA libraries**

Hands-on time: **Less than 30 minutes**

Walk-away time: **3 hours**


Throughput: **Flexible, 4 – 24 rxns in a single run**

Chemistry: **DNA Library Prep for Illumina (LP0002)**




The Volta Labs Callisto Sample Prep System offers robust and versatile library preparation through preloaded apps for high-quality short-read sequencing.

Workflow steps



The Callisto™ Sample Prep System produces high-quality Illumina libraries




Ordering information

Part number	Product	Description
C0001	Callisto™ Sample Prep System	Sample prep instrument
LP0002	DNA Library Prep for Illumina - 96 rxns	Reagents, consumables, and VoltaNodes required for 96 reactions
LP0004	DNA Library Prep for Illumina - 480 rxns, HT	Reagents, consumables, and VoltaNodes required for 480 reactions

LP0002 includes all reagents and Volta consumables required to process DNA into sequence-ready Illumina libraries.





For Research Use Only. Not for use in diagnostic procedures.

DISCLAIMER: All other product names, logos, brands, trademarks, and company names are the property of their respective owners and rights holders. Their appearance herein is not intended to and does not raise nor convey any inference, suggestion, or claim of affiliation, nor connection therewith, or endorsement thereby.