

PhotonBlade® 3 SE

advanced monopolar energy

Precision cutting. Smoke evacuation. Illumination.
All working together as one.

It can be hard to see past the smoke. Every year, 500,000 healthcare workers are exposed to harmful surgical plume.¹ That's why we're adding an advanced monopolar energy device with integrated smoke evacuation to our comprehensive portfolio.

PhotonBlade 3 SE gives you what you need all in one innovative tool.

Integrated smoke evacuation

Light and pre-stretched tubing for ease of motion and tangle-free cord management.

Ergonomic

Shape, size, grip, and styling are all improved for a better overall feel.

Light switch

Provides the user with the option to illuminate the surgical site.

Illumination

Self-contained illumination eliminates need for an additional light source.

Telescopic shaft

Fully rotatable extendable shaft reaches the small surgical pocket (2-5 in/5-12 cm).

Electrode

Delivers precise energy, now with an anti-char coating.

Standard compatibility

PhotonBlade 3 SE is compatible with the Neptune waste management system, the SafeAir compact smoke evacuator and any IEC 60601- compliant electrosurgical generator with a standard 3 prong plug. This tool also provides access to a wide variety of settings and modes across a range of full-featured ESUs.

1 Canicoba ARB. Surgical smoke and occupational health. Ann Transl Med. 2022 Dec;10(24):1303. doi: 10.21037/atm-22-5631. PMID: 36660725; PMCID: PMC9843373.

Surgical Technologies

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. We do not dispense medical advice and recommend that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate Stryker's products. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: PhotonBlade and Stryker. All other trademarks are trademarks of their respective owners or holders.

ST-PB3-BROC-1589000_REV-0 | Copyright © 2025 Stryker

Stryker Instruments
1941 Stryker Way
Portage, MI 49002
[stryker.com](https://www.stryker.com)

safeOR.com