


GINGISTAR

to Prevent Peri-implantitis



Specification

Product Name	GINGISTAR
Dimension / Weight	D 43 x W 30.2 x H 282 mm / 112g
Speed of Mode (RPM)	Weak Mode : 15,000 rpm Strong Mode : 20,000 rpm
Charging time	3.5 Hours
Number of usage When fully charged	60 times (10 min per use)
Components	Main Body + C-Type Cable Charging Adapter + Tip (10ea) + Tube (10ea) + 10cc Syringe (5ea)

GINGISTAR

Brush to prevent peri-implantitis



Product features

- Remove plaque or inflammation formed around the gum
- Prevention of Peri-implantitis
 - * Use when implanted patient(s) is recalled, to cleanse implant
- Longer battery life and improved water leakage



Inflammation occurs around abutment in the mandibular molar gum



Disinfection liquid being injected using GINGISTAR with vibration



Treating inflammation around the abutment using GINGISTAR



Inflammation is healed after using GINGISTAR

GINGISTAR Tip



2S Type



3S Type

Model name	GINGISTAR Tip (2S)	GINGISTAR Tip (3S)
Tip Color	White	Blue
Diameter	0.8 mm	0.6 mm
External Diameter	Front : 2.0mm End : 3.0mm	2.0 mm
Classification	2S Type	3S Type

- 1 Box = Tip (50ea) + Tube (50ea)