



IS-III active S-narrow




Super Narrow for Anterior Mandible (Narrow Ridge)



Diameter	Length	Product name
Ø3.2	8.5mm	IS3SN3008AP
	10.0mm	IS3SN3010AP
	11.5mm	IS3SN3011AP
	13.0mm	IS3SN3013AP


* Cover screw not included.

IS-III active S-Narrow Features




S.L.A. Surface

The new S.L.A Surface with
40% greater sur-face area and
50% more cell adhesion
promotes osseointegration



Straight Body

Easy to adjust
implantation depth



Deep & Wide Pitch

Optimum Pitch for
osseointegration
(0.9Pitch)




Wide Cutting Edge

Wide cutting edge and
enlarged surface area
enhances initial fixation and
offers clinicians more stable
implant placement




• M2.0 Connection

- One Piece Type
- Cover Screw
- Healing Abutment



• Internal Conical Seal 11° & 2/1 Hex



• M1.6 Connection

- Two Piece Type
- PickCap Impression Coping
- Temporary Abutment
- Cement Abutment
- Angled Abutment



S-narrow Prosthesis

- IS-III active S-narrow
- Lab Analog
- Cover Screw
- Healing Abutment
- PickCap Impression Coping
- Temporary Abutment
- Cement Abutment
- Angled Abutment

