

# iStatist

# HEPATITIS B

## Surface Antigen Test




## INTRODUCING

# the iStatist Hepatitis B Surface Antigen Test


- ✓ Accurate and easy-to-read results
- ✓ Results in just 15 minutes
- ✓ Simple-to-use, easy-to-store, and portable

# TEST PROCEDURE


## TEST SET UP




## BLOOD SAMPLE COLLECTION



## TESTING OF SAMPLE



## READ TEST RESULTS



## PRODUCT INFORMATION

INFORMATION TYPE	PRODUCT DETAILS
Method	Lateral Flow
Time To Results	15 Minutes
Storage Conditions	2–30 °C (35.6–86 °F)
Shelf Life	12 Months <sup>1</sup>
Sample Type	Whole Blood, Fingerstick Whole Blood, Serum or Plasma

Notes: <sup>1</sup>Shelf-life to be extended.

## ORDERING INFORMATION

PRODUCTS	
901136140030-50	iStatis HBsAg, Singles, CA, 50 (Individual tests)
901136140031-100	iStatis HBsAg, Singles, CA, 100 (Individual tests)
901149140032-50	iStatis HBsAg, Bulk, CA, 50 (Bulk tests)
901149140033-100	iStatis HBsAg, Bulk, CA, 100 (Bulk tests)

## ACCURACY

SAMPLE TYPE	PPA	NPA
Fingerstick Capillary Blood	100.0% [99.1% – 100.0%]	100.0% [99.6% - 100.0%]
Venous Whole Blood	99.5% [98.2% – 99.9%]	100.0% [99.6% - 100.0%]
Plasma	99.5% [98.2% – 99.9%]	100.0% [99.6% - 100.0%]
Serum	99.5% [98.2% – 99.9%]	100.0% [99.6% - 100.0%]

- ✓ Built-In Control
- ✓ All Support Materials Included
- ✓ Easy-To-Use and Clear Results
- ✓ High Accuracy and Reliability