

SpaceFile® Interproximal Reduction Kit

- Precision Swiss-made, electro-chemical bonded diamond files accomplish a wide range of tasks including breaking contact, reduction, contouring, fine trimming, and finishing
- Double-sided or single-sided (left or right)
- Can be used manually two ways: hand-held and with removable grip
- Files, grip and inner tray are autoclavable



Assorted Packs

Includes 8 Files (1 of each thickness), 1 Serrated Contact Point Saw, and 1 Removable Grip

DOUBLE-SIDED	SINGLE-SIDED LEFT	SINGLE-SIDED RIGHT
40034	40034-LS	40034-RS

Single Packs

Includes 8 files (all same thickness) and 1 Removable Grip

		APPROXIMATE THICKNESS		
Polishing A final step in tooth surface restoration.		40035-G	40035-LS-G	40035-RS-G
		0.13 mm	0.11 mm	0.11 mm
Finishing For finishing prior to final polish.		40035-P	40035-LS-P	40035-RS-P
		0.16 mm	0.13 mm	0.13 mm
XX Fine Pre-polish and finish reduced areas. Also use as initial strip for excessive crowding.		40035-W	40035-LS-W	40035-RS-W
		0.18 mm	0.16 mm	0.16 mm
X Fine Pre-polish stripped areas. Also use as initial strip for crowded cases.		40035-Y	40035-LS-Y	40035-RS-Y
		0.21 mm	0.19 mm	0.19 mm
Fine Use to initiate reduction in crowded cases. Also use to contour after reduction.		40035-R	40035-LS-R	40035-RS-R
		0.29 mm	0.22 mm	0.22 mm
Medium Use for reduction with some crowding or when some space exists. Smaller grain can achieve results before using Coarse strip.		40035-BL	40035-LS-BL	40035-RS-BL
		0.37 mm	0.27 mm	0.27 mm
Coarse Primary working strip for IPR. Also use to trim and contour.		40035-GR	40035-LS-GR	40035-RS-GR
		0.46 mm	0.33 mm	0.33 mm
X Coarse Use for maximum reduction at start of IPR.		40035-BLK	40035-LS-BLK	40035-RS-BLK
		0.54 mm	0.41 mm	0.41 mm
Serrated Contact Point Saw Use to break contact points on enamel, composite, amalgam, porcelain or gold.		40035-CS 0.34 mm		

SpaceFile® Removable Grip

40035-GRIP 3-Pack



Colors are for illustration purposes only. Actual color of product may vary.