



GCX Channel

for owner provided GCX style swing arms, integrates the GCX Channel to a single equipment rail, for use with installed Modular Channel. GCX hardware factory attached to rail mounting adapter. Does not include patient monitor and swing arm.

Item Number:
GCX149-MODULAR
Requires user defined
JSN number

Dimensions:
5.5"H x 2.75"W x 1.5"D
Weight: 2lbs

Item Description:

- For use with installed Hillrom Vertical Channel on Headwalls / Power Columns
- Standard adapters are plastic with embedded anti-microbial compound.
- Adapter Plates are aluminum with an anti-microbial powder coat painted finish (AMF).
- Some accessories which require heavy-duty capacity will be provided with extruded aluminum adapter(s) and/or wall stand offs for added strength.
- Some plates have multiple holes to attach future products from different vendors without the need to purchase new plates each time.
- Some plates allowing adaptability to user demand, new codes or new vendor products without going to the wall and adjusting the rail itself.
- Alternate mounting kits available for other manufacturer's devices. Please contact your Sales Representative.
- Accessories now available for Modular, Hill-Rom/Wittrock and AMICO Style Rails.