



RODIN®
3D Resin

A Pac-Dent Brand

Glaze N₂-Free

Light-Curable Naturalizing Resin

Glaze N₂-Free is a universal, light-curable sealant that neutralizes the oxygen inhibition layer, creating a smooth, enamel-like finish with less post-processing time.


BENEFITS & FEATURES

- **Superior Esthetics**
Smooth, glossy, enamel-like finish
- **Seal and Repair**
Seals surfaces and corrects minor imperfections
- **Versatile Application**
Compatible with multiple materials and intraoral use




Light-curable formula
for fast, efficient results.

Simulated Clinical Longevity (350,000 cycles - approx. 2.5 years)



Gloss Before Aging



*Data provided by Dr. Giordano of Boston University

MULTIPLE APPLICATIONS



Printed Materials



Zirconia



PMMA



Hybrid Ceramics

ORDER INFORMATION:

24031 Glaze N₂-Free, 5 mL Bottle

24015 Glaze N₂-Free, 15 mL Bottle



Pac-Dent, Inc.
670 Endeavor Circle, Brea, CA 92821
www.rodin-3d.com | www.pac-dent.com

Customer Support
(909) 839-0888
info@pac-dent.com

Technical Support
(714) 990-5488
rodin@pac-dent.com



Made in USA