

RADIESSE DILUTION REFERENCE CHART

DILUTION TYPE	RADIESSE (CAHA): SALINE/LIDOCAINE RATIO	TOTAL VOLUME PER 1.5 ML SYRINGE	USE AREAS
STANDARD (UNDILUTED)	1.5 mL Radiesse+	1.5 mL	DEEP FACIAL VOLUMIZATION (CHEEKS, JAWLINE, CHIN) LAYERS 4 OR 5 FACE
DILUTE (1:1)	1.5 mL Radiesse Classic + 1.0 mL saline/ 0.5 mL lidocaine 2%	3.0 mL	FACE (SUBDERMAL BIO-STIMULATION), LATERAL FACE, PREAURICULAR
HYPERDILUTE (1:2)	1.5 mL Radiesse Classic + 2.5 mL BNS /0.5 mL lidocaine 2%	4.5 mL	NECK, DECOLLETAGE, LOWER FACE, SKIN TIGHTENING
HYPERDILUTE (1:3)	1.5 mL Radiesse + 4.0 mL BNS + 0.5 mL lidocaine 2%	6.0 mL	CHEST (DECOLLETE), NECK, FACE OR BODY AREAS WITH THIN SKIN OR ADVANCED LAXITY

INJECTION TECHNIQUE TIPS BY AREA

FACE (BIOSTIMULATOR)	CHEST (DÉCOLLETÉ)	NECK
DILUTE 1:1 (IDEAL NEOCOLLAGENSIS) 1:2 THINNER DERMIS 1:3 SKIN QUALITY	HYPERDILUTE 1:3	HYPER-DILUTE 1:2 TO 1:3 DILUTION VARIES BASED ON DERMAL THICKNESS
SUBDERMAL FANNING OR GRID PATTERN	VERY SUPERFICIAL DERMAL PLANE	SUPERFICIAL DERMAL/SUBDERMAL PLANE
IMPROVES TEXTURE, MILD TIGHTENING, AND FINE LINES	TARGETS CREPINNESS AND FINE LINES FROM SUN DAMAGE	USE CANNULA WITH LINEAR THREADING/FANNING



LAYER 4
STANDARD
DILUTION

LAYER 2
1:1

UNDILUTE
RADIESSE

RADIESSE
MAY ALSO HYBRID
1:1
1:2
1:3

RADIESSE
1:2 OR 1:3

HYBRID
RADIESSE
BELOTERO
PLUS 0.5 NS

RADIESSE
1:3