

Reign Medical Clench™ Compression Staple



SPARTAN®
M E D I C A L

Single-Use Sterile Nitinol Compression Staple System




07x07




09x07




11x10




15x12



18x18




20x20



Technical Features


- No heat activation required
- Constant compression
- Smooth profile
- All-in-one instruments & implant kit
- Sterile, single use kits
- Reloadable instrument design
- Multiple tamp options
- Robust instrumentation



SIZE	WIDTH (W) mm	LENGTH (L) mm	THICKNESS mm	COMPRESSION mm
7mm x 7mm	7mm	7mm	1.5	1.9
9mm x 7mm	9mm	7mm		1.9
11mm x 10mm	11mm	10mm		2.8
15mm x 12mm	15mm	12mm	2.3	3.4
18mm x 18mm	18mm	18mm		5.3
20mm x 20mm	20mm	20mm		5.9



Surgical Procedure



VA Ref. No. DoD Ref. No. Description

CLIN0645

DLIN3140

Staple, Sterile Kit (All-inclusive)

LEARN MORE:

Contact Spartan Medical Inc.



(888) 240-8091



CS@SPARTANMEDICAL.COM



SPARTAN MEDICAL INC.

1445 Research Blvd, Suite 520,
Rockville, MD 20850

spartanmedical.com

UEI: LCJBD88KQ8G3

DUNS: 807438341 | CAGE:65VT5

