

Malleolar 4.0 Cannulated Screw Surgical Procedure Kit



SPARTAN®
M E D I C A L

Low Profile

Reverse-cutting flutes

Spherical underside

Self-tapping

1/3 Partial Thread



T15 Hexalobe



All implants made from
medical-grade titanium alloy

MALLEOLAR 4.0

CANNULATED SCREW KITS

38MM CANNULATED SCREW KITS

VA Ref. No. / DoD Ref. No.

Procedure Kit

CLIN01294 / DLIN4008

Description	Qty
1. 4.0 Cannulated Screws (38mm)	2
2. 7.0mm Washer	2
3. 1.25mm x 150mm Smooth Guide Wire	2
4. 1.25mm x 150mm Threaded Guide Wire	2
5. 2.7mm Cannulated Drill Bit	1
6. Parallel Wire Guide	1
7. Direct Measurement Device	1
8. Cannulated T15 Screwdriver	1
9. Canulated T15 Driver Shaft, AO	1

40MM CANNULATED SCREW KITS

VA Ref. No. / DoD Ref. No.

Procedure Kit

CLIN01295 / DLIN4009

Description	Qty
1. 4.0 Cannulated Screws (40mm)	2
2. 7.0mm Washer	2
3. 1.25mm x 150mm Smooth Guide Wire	2
4. 1.25mm x 150mm Threaded Guide Wire	2
5. 2.7mm Cannulated Drill Bit	1
6. Parallel Wire Guide	1
7. Direct Measurement Device	1
8. Cannulated T15 Screwdriver	1
9. Canulated T15 Driver Shaft, AO	1



All items are available for procurement through National VA Surgical Implants Next Generation (SING) IDIQ Contract #: 36C10G23D0066, Defense Logistics Agency (DLA) Troop Support's Distribution and Pricing Agreement (DAPA) Contract #: SP0200-19-H-0051, or through the Medical/Surgical Prime Vendors (MSPVs) - Cardinal Health PVON 118845 and Owens & Minor PVON 5967. **Electronic Data Interchange (EDI) orders and Government Purchasing Cards (GPCs) are accepted.**

Please visit our website at www.spartanmedical.com for additional product information.

SINGLE-USE.
STERILE PACKED.
READY TO USE.

Learn More:



[Link](#)

LEARN MORE:

Contact Spartan Medical Inc.

(888) 240-8091

CS@SPARTANMEDICAL.COM



SPARTAN MEDICAL INC.

1445 Research Blvd, Suite 520,

Rockville, MD 20850

spartanmedical.com

UEI: LCJBD88KQ8G3

DUNS: 807438341 | CAGE:65VT5

