

Straight Compression Plating System Surgical Procedure Kit



SPARTAN®
M E D I C A L

CORTICAL SCREWS



LOCKING SCREWS



Rounded plates and smooth edges
lessen soft tissue irritation

Protruding boss elevates
plate for reduced contact

K-wire holes for provisional fixation

Slots enable plate
positioning and compression

Threaded holes designed to
accommodate locking screws

Low profile

All implants made from
medical-grade titanium alloy

3.5MM, 6-HOLE KITS

VA Ref. No. / DoD Ref. No.

Procedure Kit	CLIN01296 / DLIN4010
Description	Qty
1. Straight Compression Plate (6-hole)	1
2. 2.8/3.5mm Double-Ended Drill Guide	1
3. 2.8mm Drill Bit	1
4. Small Depth Gauge	1
5. Solid T15 Screwdriver	1
6. 3.5mm Cortical Screws (12mm, 14mm, 16mm, 18mm, 22mm)	2 (each)

3.5MM, 8-HOLE KITS

VA Ref. No / DoD Ref. No.

Procedure Kit	CLIN01297 / DLIN4011
Description	Qty
1. Straight Compression Plate (8-hole)	1
2. 2.8/3.5mm Double-Ended Drill Guide	1
3. 2.8mm Drill Bit	1
4. Small Depth Gauge	1
5. Solid T15 Screwdriver	1
6. 3.5mm Cortical Screws (12mm, 14mm, 16mm, 18mm, 22mm)	2 (each)



ANCILLARY SCREWS*

3.5mm Cortical

8-40mm - 2mm Increments

40-50mm - 5mm Increments

3.5mm Locking

10-40mm - 2mm Increments

40-50mm - 5mm Increments

**Individually packaged as separate ancillary items*

All items are available for procurement through National VA Surgical Implants Next Generation (SING) IDIQ Contract #: 36C10G23D0066, Defense Logistics Agency (DLA) Troop Support's Distribution and Pricing Agreement (DAPA) Contract #: SP0200-19-H-0051, or through the Medical/Surgical Prime Vendors (MSPVs) - Cardinal Health PVON 118845 and Owens & Minor PVON 5967. **Electronic Data Interchange (EDI) orders and Government Purchasing Cards (GPCs) are accepted.**

Please visit our website at www.spartanmedical.com for additional product information.

**SINGLE-USE.
STERILE PACKED.
READY TO USE.**

Learn More:



[Link](#)

LEARN MORE:

Contact Spartan Medical Inc.

(888) 240-8091

CS@SPARTANMEDICAL.COM



SPARTAN MEDICAL INC.

1445 Research Blvd, Suite 520,
Rockville, MD 20850

spartanmedical.com

UEI: LCJBD88KQ8G3

DUNS: 807438341 | CAGE:65VT5

