

Straight Fibula™ Plating System Surgical Procedure Kit



SPARTAN®
M E D I C A L

CORTICAL SCREWS



LOCKING SCREWS




Rounded plates and smooth edges lessen soft tissue irritation

Slots enable plate positioning and compression

Threaded holes designed to accommodate locking screws

Low profile



All implants made from medical-grade titanium alloy

3.5MM PLATING SYSTEM KIT

VA Ref. No. / DoD Ref. No.

Procedure Kit

CLIN01298 / DLIN4012

Description	Qty
1. Straight Fibula Plate (6-hole)	1
2. 2.8/3.5mm Double-Ended Drill Guide	1
3. 2.8mm Drill Bit	1
4. 3.5mm Drill Bit	1
5. Small Depth Gauge	1
6. Solid T15 Screwdriver	1
7. 3.5mm Cortical Screws (12mm, 14mm, 16mm, 18mm, 22mm)	2 (each)



ANCILLARY SCREWS*

3.5mm Cortical

8-40mm - 2mm Increments 40-50mm

- 5mm Increments

3.5mm Locking

10-40mm - 2mm Increments

40-50mm - 5mm Increments

*Individually packaged as separate ancillary items

All items are available for procurement through National VA Surgical Implants Next Generation (SING) IDIQ Contract #: 36C10G23D0066, Defense Logistics Agency (DLA) Troop Support's Distribution and Pricing Agreement (DAPA) Contract #: SP0200-19-H-0051, or through the Medical/Surgical Prime Vendors (MSPVs) - Cardinal Health PVON 118845 and Owens & Minor PVON 5967. **Electronic Data Interchange (EDI) orders and Government Purchasing Cards (GPCs) are accepted.**

Please visit our website at www.spartanmedical.com for additional product information.

**SINGLE-USE.
STERILE PACKED.
READY TO USE.**

Learn More:



[Link](#)

LEARN MORE:

Contact Spartan Medical Inc.

(888) 240-8091

CS@SPARTANMEDICAL.COM



SPARTAN MEDICAL INC.

1445 Research Blvd, Suite 520,
Rockville, MD 20850

spartanmedical.com

UEI: LCJBD88KQ8G3

DUNS: 807438341 | CAGE:65VT5

