

SYNTHETIC BONE GRAFT - DENTAL




SPARTAN[®]
M E D I C A L

in|road[®]

DENTAL SYNTHETIC BONE GRAFT

**A patient-first,
surgeon-friendly
bone graft**

InRoad[®] emulates the biological microenvironment of bone without donor risk and moderates elaborate clinical steps without complication.
It's the bone graft of the future.




Features & Benefits

- 100% synthetic; no risk of disease transmission or donor site morbidity
- No immunogenic concerns
- Interconnected multi-porous structure promotes wettability
- Trabecular-like structure facilitates cell migration and bone regeneration
- Easy to fill bony defects by the nature of granule type bone graft




OSTEOGENE[™]


Interconnected Multi-porous Structure



■ Primary Pores



■ Secondary Channel




■ Tertiary Holes

InRoad® Dental Synthetic Bone Graft & Common Uses

OEM #	IDIQ #	VOLUME	GRANULE SIZE
HA90D05S	CLIN01207	0.5 cc	*0.3 ~ 1.0 mm
HA90D10S	CLIN01208	1.0 cc	
HA90D20S	CLIN01209	2.0 cc	
HA90D30S	CLIN01210	3.0 cc	
HA90D05L	CLIN01211	0.5 cc	0.8 ~ 1.8 mm
HA90D10L	CLIN01212	1.0 cc	
HA90D20L	CLIN01213	2.0 cc	
HA90D30L	CLIN01214	3.0 cc	



Extraction sites



Alveolar bony defects



Sinus floor elevations



Other bone defects

NOTE: Composition: HA 90 ± 5% / β-TCP 10 ± 5%

Available for procurement through VA IDIQ Contract: 36C10G23D0066 and *DoD ECAT Contract #: SPE2DF-26-D-0008



© 2022 Osteogene Tech Corp. All rights reserved.

InRoad® represents an FDA-cleared (K202675) synthetic dental bone graft intended for bone defect surgical procedures as outlined in the Instructions for Use. InRoad® and other trademarks and logos associated with InRoad® are exclusively owned by Osteogene Tech Corp and its related and affiliated companies. The use or publishing of the mark InRoad® is limited to the licensees of Osteogene Tech Corp, which is the sole and exclusive manufacturer and supplier of the InRoad® products. The InRoad® products must be used by licensed professionals (Rx Only) and sold or distributed either directly to them or through authorized distributors of Osteogene Tech Corp. Osteogene Tech Corp is not responsible or liable for any outcome resulting from any failure to strictly comply with the Instructions for Use. Results may vary. More information is available at www.osteogene.com.

IR001-REV-C 04/26

LEARN MORE:

Contact Spartan Medical Inc.

(888) 240-8091

CS@SPARTANMEDICAL.COM



SPARTAN MEDICAL INC.

7537 Lindbergh Drive
Gaithersburg, MD 20789
spartanmedical.com

DUNS: 807438341 | CAGE:65VT5

