

AlloFuse®


DEMINERALIZED BONE MATRIX

You have options,
select what's right
for your patient.

AlloFuse Demineralized Bone Matrix (DBM) putties, paste and gel are bone graft substitutes providing exceptional handling characteristics in moldable and flowable formats.



AlloFuse® Putty
DEMINERALIZED BONE MATRIX



AlloFuse® Plus Putty
DEMINERALIZED BONE MATRIX
WITH CANCELLOUS



AlloFuse® Plus Paste
DEMINERALIZED BONE MATRIX
WITH CANCELLOUS



AlloFuse® Gel
DEMINERALIZED BONE MATRIX

THE BENEFITS OF ALLOFUSE DBM PUTTY, PASTE AND GEL

Osteoinductivity is Tested on Every Lot

Ensures osteoinductive potential

Stored at Ambient Temperature

Simplifies storage and handling

Demineralized Bone with Carrier

Provides exceptional handling characteristics

Moldable and Flowable Options

Maintains desired shapes and simplifies placement of graft

Available with Added Ground Cancellous Bone

Delivers additional osteoconductive potential

Sterility Assurance Level of 10⁻⁶

Minimizes risk of infection

AlloFuse® Putty

DEMINERALIZED BONE MATRIX

VOLUME	CONDITION	STORAGE	OEM #	IDIQ REF NO	DAPA REF NO
1 cc	Sterile	Ambient	90038001	CLIN0913	DLIN3635
2.5 cc	Sterile	Ambient	90038002	CLIN0914	DLIN3636
5 cc	Sterile	Ambient	90038005	CLIN0915	DLIN3637
10 cc	Sterile	Ambient	90038010	CLIN0916	DLIN3638

Packaged in a jar

AlloFuse® Plus Putty

DEMINERALIZED BONE MATRIX WITH CANCELLOUS

VOLUME	CONDITION	STORAGE	OEM #	IDIQ REF NO	DAPA REF NO
5 cc	Sterile	Ambient	90338005	CLIN0923	DLIN3819
10 cc	Sterile	Ambient	90338010	CLIN0924	DLIN3646

Packaged in a jar

AlloFuse® Plus Paste

DEMINERALIZED BONE MATRIX WITH CANCELLOUS

VOLUME	CONDITION	STORAGE	OEM #	IDIQ REF NO	DAPA REF NO
1 cc	Sterile	Ambient	90238001	CLIN0920	DLIN3818
3 cc	Sterile	Ambient	90238003	CLIN0921	DLIN3643
8 cc	Sterile	Ambient	90238008	CLIN0922	DLIN3644

Packaged in a syringe

AlloFuse® Gel

DEMINERALIZED BONE MATRIX

VOLUME	CONDITION	STORAGE	OEM #	IDIQ REF NO	DAPA REF NO
1 cc	Sterile	Ambient	90138001	CLIN0917	DLIN3639
5 cc	Sterile	Ambient	90138005	CLIN0918	DLIN3640
10 cc	Sterile	Ambient	90138010	CLIN0919	DLIN3641

Packaged in a syringe

- Coetzee J.C. Use of demineralized bone matrix gel in foot surgery: A prospective randomized controlled study comparing demineralized bone matrix to autologous bone graft. *AlloSource White Paper*. 2012; M8S0040.04.
- Walsh W.R., Oliver R, Yu Y, Bell D.J., Pak P, Russel N. Demineralized bone matrix provides equivalent results to autograft in standard posterolateral fusion model in adult rabbits. *AlloSource White Paper*. 2012; M8S0039.00.

AlloFuse®

PERFORMANCE AND HANDLING

Trusted and used in over 400,000 procedures since 2006

Shown to support bone healing in clinical and pre-clinical studies^{1,2}

Exceptional handling and placement characteristics

- Easily moldable at cooler temperatures (like the OR)
- Becomes thicker at warm temperatures (like the human body) to stay in place
- Resists irrigation

CLIN numbers are available on VA IDIQ Contract: 36C10G23D0066

DLIN numbers are available on DLA DAPA Contract: SP0200-19-H-0051 and Prime Vendors Cardinal Health (PVON 118845) and Owens & Minor (PVON 5967)

LEARN MORE:

Contact Spartan Medical Inc.



(888) 240-8091



CS@SPARTANMEDICAL.COM



SPARTAN MEDICAL INC.

7537 Lindbergh Drive
Gaithersburg, MD 20879
spartanmedical.com
UEI: LCJBD88KQ8G3
DUNS: 807438341 | CAGE:65VT5



Distributed by Spartan Medical. Manufactured by AlloSource.

AlloFuse Demineralized Bone Matrix putties, paste and gel are a 510k cleared medical device. AlloSource is registered with the FDA as a tissue establishment and accredited by the Association for Advancing Tissue and Biologics.

© Spartan Medical Inc. 2026 | AlloFuse® is a registered trademark of AlloSource

