

AlloFuse[®]

FIBERS AND FIBER BOATS^{*}


DEMINERALIZED BONE

You have options, select what's right for your patient.


AlloFuse Fibers and Fiber Boats can be used in spinal, foot and ankle or other orthopedic procedures to fill boney voids. They provide an osteoconductive scaffold with osteoinductive potential, and may be used with biologics or autograft tissue as directed by the surgeon.



AlloFuse[®] FIBERS
DEMINERALIZED BONE



AlloFuse[®] MICRO FIBERS
DEMINERALIZED BONE



AlloFuse[®] FIBER BOAT
DEMINERALIZED BONE

THE BENEFITS

OF ALLOFUSE FIBERS AND FIBER BOATS

- ⊕ **MADE OF 100% DEMINERALIZED CORTICAL BONE**
Contains no carrier
- ⊕ **MADE OF 100% DEMINERALIZED CORTICAL BONE**
Contains no carrier
- ⊕ **MAY BE USED WITH BIOLOGICS OR AUTOGRAFT TISSUE**
Developed to support bone growth
- ⊕ **MOLDABLE AND FORMABLE**
Users control consistency to maintain desired shapes
- ⊕ **STERILITY ASSURANCE LEVEL OF 10⁻⁶**
Minimizes risk of infection
- ⊕ **OSTEOINDUCTIVITY IS TESTED ON EVERY LOT**
Ensures osteoinductive potential
- ⊕ **STORED AT AMBIENT TEMPERATURE**
Simplifies storage and handling

AlloFuse®

FIBERS AND FIBER BOATS*

DEMINERALIZED BONE

YOUR EXPERTISE, YOUR DECISION.

ALLOFUSE FIBERS

STORAGE: AMBIENT

VOLUME	OEM #
1 cc	88231001
5 cc	88231005
10 cc	88231010



Fiber Strands

- Form into compact shapes to fill **large boney voids**.
- Optimally designed for procedures requiring higher volumes of allograft.

ALLOFUSE MICRO FIBERS

STORAGE: AMBIENT

VOLUME	OEM #	IDIQ REF NO	DAPA REF NO
1 cc	88241001	CLIN01168	DLIN3950
5 cc	88241005	CLIN01169	DLIN3951
10 cc	88241010	CLIN01170	DLIN3952



Fiber Strands

- Form into tightly compacted shapes to fill **small boney voids**.
- Eliminate cutting or mincing, to save time and reduce unnecessary work.

ALLOFUSE FIBER BOAT

STORAGE: AMBIENT

SIZE	LENGTH	WIDTH	HEIGHT	VOLUME	OEM #
Small	50 mm	25 mm	6 mm	5 cc	88031050
Large	100 mm	25 mm	6 mm	10 cc	88031100



Fiber Boat

- May be used with biologics or autograft tissue to replace patient bone as directed by the surgeon.
- Small and large sizes available to fill voids or span across bone structures.

CLIN numbers are available on VA IDIQ Contract: 36C10G23D0066

DLIN numbers are available on DLA DAPA Contract: SP0200-19-H-0051 and Prime Vendors Cardinal Health (PVON 118845) and Owens & Minor (PVON 5967)

LEARN MORE:

Contact Spartan Medical Inc.



(888) 240-8091



CS@SPARTANMEDICAL.COM



SPARTAN MEDICAL INC.

7537 Lindbergh Drive
Gaithersburg, MD 20879
spartanmedical.com
UEI: LCJBD88KQ8G3
DUNS: 807438341 | CAGE:65VT5



Distributed by Spartan Medical. Manufactured by AlloSource.

AlloFuse® Fibers and Fiber Boats are regulated by the FDA under 21 CFR Part 1271 Human Cells, Tissues, and Cellular and Tissue-Based Products (HCT/PS). AlloSource is registered with the FDA as a tissue establishment and accredited by the Association for Advancing Tissue and Biologics.

© Spartan Medical Inc. 2026 | AlloFuse® is a registered trademark of AlloSource

