

ANTHEM® MADE IN AMERICA FOR EVERYDAY AMERICAN HEROES



Anthem®

Aurelia® Anthem® Gloves are changing the game. **Introducing the first-ever American-made glove from Aurelia®**, built with pride, precision, and performance. From healthcare to dental to sanitation, Anthem® delivers what professionals demand most.

Available Through ECAT Contract #: SPE2DF-26-D-0008

FEATURES

4.5mil Finger-Textured Construction
Powder-Free Nitrile Exam Gloves
Chemo and Fentanyl Rated
Beaded Cuff
Amidextrous
AQL 1.5

INSTRUCTION FOR USE

NON-STERILE. FOR SINGLE USE ONLY. Before putting on gloves hands must be clean and dry. Gloves do not replace the need for hand hygiene. Hands must be washed every time gloves are removed to avoid contamination. Gloves are for single use. Change gloves if moving from a contaminated area to a non-contaminated area. Be sure to wear gloves that fit properly, gloves should fit snugly around wrists and hands. Different sizes are available.

STORAGE

Store in a cool, dry place. Avoid excessive heat (40°C, 104°F). Shield open boxes from direct sun and fluorescent lighting. To order, please refer to a specific size and mention the product code.

TENSILE TEST RESULTS (unaged)

	Tensile Strength (MPa)	Elongation at break (%)
XS	30.0	680
S	31.1	680
M	31.0	680
L	31.2	680
XL	31.3	680
XXL	31.3	680
ASTM D6319	14 MPa minimum	500% minimum

AGED TENSILE TEST RESULTS (7 Days at 70° C)

	Tensile Strength (MPa)	Elongation at break (%)
XS	30.0	620
S	30.0	620
M	30.0	620
L	30.0	620
XL	30.0	620
XXL	30.0	620
ASTM D6319	14 MPa minimum	620% minimum

DESCRIPTION

Pallet	OEM #	Size	Qty	Length	Thickness	
					Palm (mils/mm)*	Finger (mils/mm)*
17765-P	17765	XS	10 boxes of 100 gloves/case	9.4 / 240 8.7 / 220 min.	4.0 / 0.10 1.97 / 0.05 min.	4.5 / 0.11 1.97 / 0.05 min
17766-P	17766	S	10 boxes of 100 gloves/case	9.4 / 240 8.7 / 220 min.	4.0 / 0.10 1.97 / 0.05 min.	4.5 / 0.11 1.97 / 0.05 min
17767-P	17767	M	10 boxes of 100 gloves/case	9.4 / 240 8.7 / 220 min.	4.0 / 0.10 1.97 / 0.05 min.	4.5 / 0.11 1.97 / 0.05 min
17768-P	17768	L	10 boxes of 100 gloves/case	9.4 / 240 8.7 / 220 min.	4.0 / 0.10 1.97 / 0.05 min.	4.5 / 0.11 1.97 / 0.05 min
17769-P	17769	XL	10 boxes of 90 gloves/case	9.4 / 240 8.7 / 220 min.	4.0 / 0.10 1.97 / 0.05 min.	4.5 / 0.11 1.97 / 0.05 min
17760-P	17760	XXL	10 boxes of 80 gloves/case	9.4 / 240 8.7 / 220 min.	4.0 / 0.10 1.97 / 0.05 min.	4.5 / 0.11 1.97 / 0.05 min

*Thickness (single-wall)

NOTE: Shipping included on all pallet orders (70 cases = 1 pallet)



LEARN MORE:

Contact Spartan Medical Inc.

(888) 240-8091

CS@SPARTANMEDICAL.COM



SPARTAN MEDICAL INC.

7537 Lindbergh Drive
 Gaithersburg, MD 20879

spartanmedical.com

UEI: LCJBD88KQ8G3

DUNS: 807438341 | CAGE:65VT5





SPARTAN®
M E D I C A L

AURELIA® GLOVES

Available Through ECAT Contract #: SPE2DF-26-D-0008

The **AURELIA Anthem Gloves**, a Supermax Healthcare Inc. brand, is proudly manufactured in the U.S. by Maxter Healthcare Inc., a vertically integrated manufacturer of high-quality nitrile examination gloves located in Rosharon, Texas. Spartan Medical (SDVOSB) is the Federal Government Supplier for the AURELIA Anthem 4.5mil Nitrile, Powder-free, Chemo and Fentanyl Rated Medical Exam Gloves. Maxter Healthcare Inc. is a board member of the American Medical Manufacturing Association (AMMA), and the MAXTER® Gloves are AMMA-certified, in compliance with the recent Executive Order – [Ensuring Truthful Advertising of Products Claiming to be Made in America \(March 13, 2026\)](#).



Core Competencies

- **US-based manufacturing of high-quality, 510(k)-cleared powder-free nitrile examination gloves.**
- Compliance with ASTM D6319 (nitrile exam glove specification), ASTM D412 (tensile/elongation), ASTM D5151 (freedom-from-holes/AQL), ISO 9001, ISO 13485, ISO 14001:2015; SMETA 7.0 and NFPA 1999.
- Chemo/fentanyl resistance tested per ASTM D6978.
- AQL ≤ 1.5 (medical grade), exceeding industry standard; sampling per ISO 2859.
- Fully Automated lines: NVIDIA AI detection; Automated Packing.
- **Six-Post leaching process to establish highest quality.**
- FDA food-contact compliant (21 CFR 177.2600).
- Available specs: 4.5 mil (patriot blue, finger-textured). XS/S/M/L/XL/XXL

Made in USA

- **100% American Manufacturing** – Rosharon, TX campus positioned to be the largest nitrile-glove producer in the U.S.
- **Secure Supply Chain – Domestic production reduces tariff, port, and geopolitical risk; shorter lead times & high fill rates.**
- Sustainability & Innovation – Solar energy, water recycling, and AI-driven automation.
- Government & Defense Eligible – **Make PPE in America Act and Buy America Act compliant**
- GO TEXAN member – rooted in Texas manufacturing and jobs.

Facilities & Capacity

- Location: Rosharon, Texas (Greater Houston)
- Scalable 215-acre campus designed for up to 72 production lines (~19B glove/year at full build-out).
- Automation: Fully dual line robotic dipping with AI vision inspection. 20 million gloves/month/line
- Current Capacity: 220 million gloves/month
- Warehouse: 100,000+ sq. ft. with national JIT distribution

Scan the QR code
or click [here](#) to
learn more:



LEARN MORE:

Contact Spartan Medical Inc.

 (888) 240-8091

 CS@SPARTANMEDICAL.COM



SPARTAN MEDICAL INC.

7537 Lindbergh Drive
Gaithersburg, MD 20879

spartanmedical.com

UEI: LCJBD88KQ8G3

DUNS: 807438341 | CAGE:65VT5



Rev D 7/26