



MAXTER® Gloves

MAXTER® MADE IN AMERICA FOR EVERYDAY AMERICAN HEROES



MAXTER® Gloves

MAXTER® Gloves are changing the game. **Introducing American-made gloves from MAXTER® Healthcare**, built with pride, precision, and performance. From healthcare to dental to sanitation, MAXTER® delivers what professionals demand most.

FEATURES

6.0mil Finger-Textured Construction
Powder-Free Nitrile Exam Gloves
Chemo and Fentanyl Rated
Beaded Cuff
Ambidextrous
AQL 1.5

INSTRUCTION FOR USE

NON-STERILE. FOR SINGLE USE ONLY. Before putting on gloves hands must be clean and dry. Gloves do not replace the need for hand hygiene. Hands must be washed every time gloves are removed to avoid contamination. Gloves are for single use. Change gloves if moving from a contaminated area to a non-contaminated area. Be sure to wear gloves that fit properly, gloves should fit snugly around wrists and hands. Different sizes are available.

STORAGE

Store in a cool, dry place. Avoid excessive heat (40°(, 104°F). Shield open boxes from direct sun and fluorescent lighting. To order, please refer to a specific size and mention the product code.

Available Through VFCP – 38 U.S.C. § 8127(d)

TENSILE TEST RESULTS (unaged)

	Tensile Strength (MPa)	Elongation at break (%)
S	30.0	525
M	30.0	525
L	30.0	525
XL	30.0	525
ASTM D6319	14 MPa minimum	500% minimum

TENSILE TEST RESULTS (aged)

	Tensile Strength (MPa)	Elongation at break (%)
S	28.0	450
M	28.0	450
L	28.0	450
XL	28.0	450
ASTM D6319	14 MPa minimum	400% minimum

DESCRIPTION

OEM #	Size	Qty	Length	Thickness	
				Palm*	Finger*
M73996	S	100/box	9.40 inches 9.05 inches min.	4.00 mil 1.97 mil min.	6.00 mil 1.97 mil min.
M73997	M	100/box	9.40 inches 9.05 inches min.	4.00 mil 1.97 mil min.	6.00 mil 1.97 mil min.
M73998	L	100/box	9.40 inches 9.05 inches min.	4.00 mil 1.97 mil min.	6.00 mil 1.97 mil min.
M73999	XL	90/box	9.40 inches 9.05 inches min.	4.00 mil 1.97 mil min.	6.00 mil 1.97 mil min.

*Thickness (single-wall)

NOTE: Shipping included on all 2 pallet orders (70 cases = 1 pallet)



LEARN MORE:

Contact Spartan Medical Inc.

(888) 240-8091

CS@SPARTANMEDICAL.COM



SPARTAN MEDICAL INC.

7537 Lindbergh Drive
Gaithersburg, MD 20879

spartanmedical.com

UEI: LCJBD88KQ8G3

DUNS: 807438341 | CAGE:65VT5



Available Through VFCP – 38 U.S.C. § 8127(d)

MAXTER[®] Gloves are proudly manufactured in the U.S. by Maxter Healthcare Inc., a vertically integrated manufacturer of high-quality nitrile examination gloves located in Rosharon, Texas; compliant with the Make PPE in America Act and Buy America Act. Spartan Medical (SDVOSB) is an authorized Federal Government Supplier for the Maxter Black 6.0mil Nitrile, Powder-free, Chemo and Fentanyl Rated Medical Exam Gloves. Maxter Healthcare Inc. is a board member of the American Medical Manufacturing Association (AMMA), and the MAXTER[®] Gloves are AMMA-certified, in compliance with the recent Executive Order – [Ensuring Truthful Advertising of Products Claiming to be Made in America](#).



Core Competencies

- **US-based manufacturing of high-quality, 510(k)-cleared powder-free nitrile examination gloves.**
- Compliance with ASTM D6319 (nitrile exam glove specification), ASTM D412 (tensile/elongation), ASTM D5151 (freedom-from-holes/AQL), ISO 9001, ISO 13485, ISO 14001:2015; SMETA 7.0 and NFPA 1999.
- Chemo/fentanyl resistance tested per ASTM D6978.
- AQL ≤ 1.5 (medical grade), exceeding industry standard; sampling per ISO 2859.
- Fully Automated lines: NVIDIA AI detection; Automated Packing.
- **Six-Post leaching process to establish highest quality.**
- FDA food-contact compliant (21 CFR 177.2600).
- Available specs: 6.0 mil (black, finger-textured). S/M/L/XL

Made in USA

- **100% American Manufacturing** – Rosharon, TX campus positioned to be the largest nitrile-glove producer in the U.S.
- **Secure Supply Chain – Domestic production reduces tariff, port, and geopolitical risk; shorter lead times & high fill rates.**
- Sustainability & Innovation – Solar energy, water recycling, and AI-driven automation.
- Government & Defense Eligible – **Make PPE in America Act and Buy America Act compliant**
- GO TEXAN member – rooted in Texas manufacturing and jobs.

Facilities & Capacity

- Location: Rosharon, Texas (Greater Houston)
- Scalable 215-acre campus designed for up to 72 production lines (~19B glove/year at full build-out).
- Automation: Fully dual line robotic dipping with AI vision inspection. 20 million gloves/month/line
- Current Capacity: 220 million gloves/month
- Warehouse: 100,000+ sq. ft. with national JIT distribution

Scan the QR code
or click [here](#) to
learn more:



LEARN MORE:

Contact Spartan Medical Inc.

 (888) 240-8091

 CS@SPARTANMEDICAL.COM



SPARTAN MEDICAL INC.

7537 Lindbergh Drive
Gaithersburg, MD 20879

spartanmedical.com

UEI: LCJBD88KQ8G3

DUNS: 807438341 | CAGE:65VT5

