

ONLINE CONSUMER TESTS		FOR SPECTRACELL USE ONLY						REQUISITION NUMBER
Purchase Date		ACD	SST	LAV	BLUE	UNST	T-UNK	For SpectraCell Use Only
LID		T-S	T-P	OI °C	NOI °C	Sample #		
Provider/Rep		Kit Type		Verify Doc	Verify Name	Check #		
Purchaser		Tube NTM						
Ship to Address		REJ Test						
		Tube #				Docu #		
								Accession #/Date Received/ Batch

PATIENT INFORMATION			
First Name*	Last Name*	MI	
Date of Birth* (MM/DD/YYYY)	Gender*	Email Address* (for result delivery)	Phone Number*
<input type="text"/>	<input type="checkbox"/> M <input type="checkbox"/> F	<input type="text"/>	<input type="text"/>

SPECIMEN INFORMATION			
▼ TO BE COMPLETED AT TIME OF BLOOD DRAW ▼			
Date of Collection	Time*	Fasting*	Phlebotomy Vendor ID
<input type="text"/>	<input type="checkbox"/> AM <input type="checkbox"/> PM	<input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="text"/>

TEST SELECTION
Micronutrient Test (MNT)
Lipoprotein Particle Profile (LPP®)
MTHFR
Telomere
Energy Optimization
Weight Optimization
Immune Health
Mental Health
AgeNexus - Biological Age Management
CardiacNexus
Cardi-A1C
BaselineNexus

PAYMENT INFORMATION	
Order #:	
Subtotal:	
Coupon Code:	Discount Amount:
Shipping:	
Grand Total:	

A COPY OF THIS FORM MUST BE RETURNED WITH SPECIMEN.

Cancellation Policy

The following applies to purchases made through spectracell.com. In order for a customer to be eligible for a cancellation or refund, customers must provide proof of purchase.

- If you request a cancellation before your order has shipped, you are eligible for a full refund.
- If it has been less than 30 days since the order was placed online and you have not sent in your sample, you are eligible for a full refund if the test kit is returned to SpectraCell unopened and unused. The customer is responsible for return shipping. (Do not use the return label included with the kit.) If the test kit is not returned, the customer will not be eligible for a refund.
- Orders over 30 days old, where the sample has not been received, cannot be cancelled and are not eligible for a refund.

FIND A DRAWSITE

! For a list of Draw Sites to have patients have their blood sample collected, scan the QR Code below.



"Find a Drawsite" link is available for US Patients. Contact a member of our Client Services Staff for drawsites outside of these territories.

MAKE SURE TO INCLUDE REQUISITION FORM INSIDE OF PATIENT'S SAMPLE KIT WHEN SHIPPING.

Collect samples Monday-Friday, except on or the day before federal holidays.

Verify the expiration date on the test tube, and label it with the **PATIENT NAME**, **DATE OF BIRTH**, and **DATE OF COLLECTION**.

It is strongly recommended that a 21-gauge needle be used for specimen collection.

Test	Fasting	Collection Tube (Volume)	Mixed by Inverting	Centrifuge Between 30-60 mins	Refrigerate Until Shipped
Micronutrient (MNT)	N/A	ACD (2)	8 x	No	N/A
LPP Plus	Yes (9-12 hrs)	SST ¹ (2)	N/A	Yes ² (3000 rpm for 15mins)	Yes
Telomere	N/A	BLU (1)	6 x	No	N/A
MTHFR	N/A	LAV (1)	6 x	No	N/A
Immune Health	N/A	ACD (2)	8 x	No	N/A
Mental Health	N/A	ACD (2)	8 x	No	N/A
Energy Optimization	N/A	ACD (2)	8 x	No	N/A
Weight Optimization	N/A	ACD (2)	8 x	No	N/A
AgeNexus	N/A	ACD (2)	8 x	No	N/A
		BLU (1)	6 x	No	N/A
CardiacNexus	Yes (9-12 hrs)	SST ¹ (2)	N/A	Yes ² (3000 rpm for 15mins)	Yes
		LAV (1)	6 x	No	N/A
Cardi-A1c	Yes (9-12 hrs)	SST ¹ (2)	N/A	Yes ² (3000 rpm for 15mins)	Yes
		LAV (1)	6 x	No	N/A
Baseline Nexus	This kit includes the MNT, LPP, Telomere, and MTHFR tests. Please follow the individual test instructions provided above.				



ACD
ACD
Solution A



SST
Serum
Separator



LAV
EDTA



BLU
Sodium
Citrate

1. Place collection tubes directly next to (2) frozen ice bricks.

2. Do NOT pour off the serum.

SPECIMEN PACKAGING – ALL SPECIMENS SHOULD BE PACKAGED IN ONE KIT

Micronutrient (MNT) Kit Only

Refer to temperature requirements carefully as it greatly affects MNT results.

1. Freeze ice brick(s) for 24 hours before drawing MNT specimen(s) if local temperatures are above 90°F. * (Use all ice bricks included in the kit)
2. Collect the specimen according to the instructions in the chart above.
3. Insert the collection tubes inside the foam holder, then place it in the biohazard bag. Put the requisition in the outside pocket of the bio-hazard bag.
4. Place ice brick(s) in the smaller compartment and the ACD tubes in the larger compartment of the Styrofoam box. Set the Styrofoam lid securely in place and close the cardboard lid. Use tape to secure the cardboard lid.
5. One kit can be utilized for multiple Micronutrient Test (MNT) samples or patients. Kindly adhere to steps 1-3 as outlined above.

Multiple Tests in one kit: How to package a kit with multiple tubes.

1. Freeze ice brick(s) for 24 hours before drawing specimen(s) requiring refrigeration. * (Use all ice bricks included in the kit)
2. Select the appropriate collection tubes for the requested test(s) according to the instructions in the chart above.
3. All specimens should be placed in the bio-hazard bags. Put the requisition(s) in the outside pocket of one of the bio-hazard bags.
4. **SST tubes:** Place all SST tubes inside the smaller compartment directly in contact with the ice brick(s), and put all other tubes in the larger compartment.
Do not place ACD tubes directly in contact with frozen ice bricks.
5. Set the Styrofoam lid securely in place and close the cardboard lid. Use tape to secure the cardboard lid.

❖ **Domestic Shipping:** *The Priority Overnight prepaid FedEx shipping label has been affixed to the SpectraCell Kit for return shipping.*

For **FedEx**, call 1-800-GoFedEx (1-800-463-3339) to schedule a pickup.

Always leave samples with a FedEx Pick up driver or inside a FedEx facility. DO NOT USE COURIER DROP BOX.

❖ **International Shipping:** Place Air Waybill, commercial invoice and CDC notice in pouch directly affixed to the SpectraCell Kit for return shipping. For **FedEx**, call 1-800-GoFedEx (1-800-463-3339) to schedule a pickup.

Always leave samples with a FedEx Pick up driver or inside a FedEx facility. DO NOT USE COURIER DROP BOX.

All specimens **MUST BE SHIPPED OVERNIGHT** on the same day they are collected.

Call SpectraCell Labs Client Services at 800.227.5227 ext. 7350 with any questions regarding specimen requirements or shipping instructions.