


POWER OF INGENUITY

Kewaunee's Evolution is a free standing workstation system that is designed to meet the changing needs of today's Laboratory and Technical environments. The column based design offers infinite height adjustment of the worksurface and shelving; as well as other column mounted accessories. Its modular design and ability to share columns helps control cost and simplifies lab layout. Three standard column heights create platforms with and without an above counter support structure. Evolution is available prepiped and pre-wired and easily incorporates with overhead utility grids or service carriers and panels. Evolution is designed for mobile, suspended, free-standing, and overhead storage.



Rear Post	Multi Height	Multi Height	Multi Height HD
Post Shape	Rectangular	Rectangular	Rectangular
Post Shape Dimensions	9"	12"	12" HD
All Colors	✓	✓	✓
Min / Max (leveling glides)	36", 54" 84"	36", 54" 84"	36", 54" 84"
Post Material	CRS/Extrusion	CRS/Extrusion	CRS/Extrusion
Prewired and Data			
Add on Horizontal Power Bar	✓	✓	✓
Fixed Horizontal Raceway	N/A		
Vertical Duplex Above & Below	Configurable up to 8A & 4B Per Post		
Prepiped Lab Services			
Vertical (Per Post)	3	3	3
Horizontal	N/A	N/A	N/A
Customizable	✓	✓	✓
Shelving			
Above	✓	✓	✓
Under	✓	✓	✓
Shelving Adjustability	Fully	Fully	1"
Configuration	Full Length (54", 84" high only)	Full Length (54", 84" high only)	Full Length (54", 84" high only)
Details			
Work Surface Height (leveling glides)	19" - 37"+	19" - 37"+	19" - 37"+
Shipping Flat Pack / Assembled (ASM)	Flat	Flat	Flat
Accepts Suspended Cabinets	✓	✓	✓
Accepts Upper Storage Cabinets	✓, 84" high only	✓, 84" high only	✓, 84" high only
Load ratings			
Single Sided Max Load	1200 lbs	1750 lbs	2600 lbs
Double Sided Max Load	2000 lbs	2500 lbs	N/A
Single Sided Worksurface	750 lbs	1000 lbs	1100 lbs
Double Sided Worksurface	500 lbs per side	600 lbs per side	600 lbs per side
Shelves Each	250 lbs	250 lbs	500 lbs
Sizes			
Starter adder	✓	✓	✓
Widths (see catalog for increments)	24"-96"	24"-96"	24"-96"
Worksurface Depth (Single sided)	24", 30" 36"	24", 30" 36"	24", 30" 36"
Worksurface Depth (Double sided)	48"-72"	48"-72"	48"-72"
Foot Options			
Leveler	✓	✓	✓
Swivel Locking Caster	✓	✓	✓
Leveling Caster & HD	✓	✓	✓
Work Tops			
Kemresin	✓	✓	✓
Phenolic	✓	✓	✓
ESD Laminate	✓	✓	✓
Plastic Laminate	✓	✓	✓
Accessories	✓	✓	✓
Certifications			
UL Listed (Factory Wired)	61010-1	61010-1	61010-1
UL Listed Power Bar (Factory Wired)	61010-1	61010-1	61010-1
SEFA Rating / Class	10/3 or 10/8 with Casters	10/3 or 10/8 with Casters	10/3 or 10/8 with Casters

FEATURES



Accessories

- >> Fully adjustable shelves
- >> CPU holder
- >> Monitor mount
- >> Keyboard
- >> Single & double lockable drawers
- >> LED light rail
- >> Bin bar
- >> Powerstrip
- >> Foot rest

Features

- >> Infinity shelving adjustments
- >> Will allow for 33" mobile storage units
- >> Allows for up to 3 service fixtures per post
- >> Data bezel optional on horizontal raceways
- >> SEFA 10 Class 7 and Class 8
- >> Capacity: Work surface 1,000 lb, Per Shelf 40lb/sqf
- >> Ships assembled and knock-down
- >> Starter and add-on assemblies
- >> Single and double sided
- >> UL 61040-1

Evolution

COLUMN BASED WORKSTATIONS




POWER OF INGENUITY



Double sided bench
Upper cabinets
Suspended and mobile cabinets




Pass through
cabinets



Workstation
• 12" Column
• 3" Configuration 2 | Category 2

Maximum Load
2500 Pounds



POWER OF INGENUITY

