

# ENTERPRISE

---

## Free-Standing Workstations



KEWAUNEE®



# Free-Standing Workstations

- Self-supporting workstations
- Single- and double-sided
- Adjustable height worksurfaces
- Integrated shelving systems
- Pre-plumbed integrated services
- Pre-wired integrated power and data
- SEFA 10, Class 7, Category 4
- UL 61010-1 listed when pre-wired

A North Carolina based company founded in 1906, Kewaunee's history of quality products and scientific expertise positions it as a recognized global leader in the design, manufacture, and installation of laboratory, healthcare, and technical furniture products. Products include steel and wood casework, fume hoods, adaptable modular systems, movable workstations, stand-alone benches, and epoxy resin worksurfaces and sinks. Kewaunee's corporate headquarters are located in Statesville, North Carolina.




# Enterprise Workstation Information



**SEFA 10 Class 7**  
(on glides)

**SEFA 10 Class 8**  
(on casters)



## Single-sided Workstation

### 2" Post

[Page 8](#)

Shelves: on front only  
 Worksurface: on front only  
 Services: in post - plumbed with hard tubing  
 Electric/Data: in horizontal chase (one side only) and below worksurface with add-on box  
 Weight Rating: 1,600lb max, 200lb per shelf

## Double-sided Workstations

### 6" Post

[Page 13](#)

Shelves: on front and back  
 Worksurface: on front and back  
 Services: in post - plumbed with flexible tubing  
 Electric/Data: in horizontal chase (both sides) and below worksurface in post  
 Weight Rating: 1,600lb max per side, 200lb per shelf



# ENTERPRISE

## Free-Standing Workstations

Kewaunee Scientific Corporation is dedicated to manufacturing high quality products for the laboratory marketplace. We offer **Total Laboratory Solutions, Innovative Designs, and Technological Expertise** with a worldwide distribution network to provide laboratories that are truly World Class.

### Laboratory Products from Kewaunee®:

ALPHA® – Adaptable Casework System  
ALPHA Columns – Expandable Workstations  
ALPHA Overhead Carrier Systems  
BasikBench® – Simple 4-Leg Benches or Tables  
Element Bench – Efficient Benches with Plumbed Posts  
Enterprise® – Free-Standing Workstations  
Evolution® – Column Based Workstations  
Kemresin® – Epoxy Resin Counters & Sinks  
Kewaunee Matrix® – Educational Stations  
Research Collection® – Steel Casework  
Signature® Series – Wood Casework  
Supreme Air® – Fume Hoods

For our complete product offering, specifications, catalogs, brochures, case studies, and dealer list, go to [kewaunee.com](http://kewaunee.com), or contact your local sales representative directly.



#### Kewaunee Scientific Corporation

P.O. Box 1842 • Statesville, NC 28687-1842  
704-873-7202 • [kscmarketing@kewaunee.com](mailto:kscmarketing@kewaunee.com) • [kewaunee.com](http://kewaunee.com)

#### Kewaunee Labway India Pvt. Ltd.

No. CA-9A, 2nd Floor, Jigani Link Road  
Jigani, Bangalore 562 106, India  
Phone: 9180-27826725 • Fax: 9180-27826724  
Web Site: [kewaunee.in](http://kewaunee.in)

#### Kewaunee Labway Asia Pte. Ltd.

194 Pandan Loop, #06-22  
Pantech Business Hub, Singapore 128 383  
Phone: 65-6773-0288 • Fax: 65-6773-2322



KEWAUNEE®

Kewaunee, Alpha, BasikBench, Enterprise, Evolution, Kemresin, Kewaunee Matrix, Research Collection, Signature, and Supreme Air are registered trademarks of Kewaunee Scientific Corporation