

# EVOLUTION

## Column Systems



KEWAUNEE®



# EVOLUTION

## ***Column Based Workstations***

- Infinitely height adjustable work surfaces and shelving
- Starter and add-on assemblies; single and double-sided
- Laminates, edge-grain maple, phenolic resin and Kemresin work surfaces available
- UL listed, column mount power bars with power switch and built-in circuit breaker
- Power, data, and service fixtures mountable in columns with pre-wired and pre-piped options
- Storage options include mobile, suspended, freestanding, and overhead
- Meets SEFA 10, Adaptable Casework Systems, Class 7 (freestanding) and Class 8 (mobile) standards





# Typical Installations





# EVOLUTION

## Column Systems

Kewaunee Scientific Corporation is dedicated to manufacturing high quality products for the laboratory marketplace. We offer **Total Laboratory Solutions, Innovative Designs, and Technological Expertise** with a worldwide distribution network to provide laboratories that are truly World Class.

### Laboratory Products from Kewaunee®:

ALPHA® – Adaptable Casework System  
ALPHA Columns – Expandable Workstations  
ALPHA Overhead Carrier Systems  
BasikBench® – Simple 4-Leg Benches or Tables  
Element Bench – Efficient Benches with Plumbed Posts  
Enterprise® – Free-Standing Workstations  
Evolution® – Column Based Workstations  
Kemresin® – Epoxy Resin Counters & Sinks  
Kewaunee Matrix® – Educational Stations  
Research Collection® – Steel Casework  
Signature® Series – Wood Casework  
Supreme Air® – Fume Hoods

For our complete product offering, specifications, catalogs, brochures, case studies, and dealer list, go to [kewaunee.com](http://kewaunee.com), or contact your local sales representative directly.



The mark of responsible forestry



### Kewaunee Scientific Corporation

P.O. Box 1842 • Statesville, NC 28687-1842  
704-873-7202 • [kscmarketing@kewaunee.com](mailto:kscmarketing@kewaunee.com) • [kewaunee.com](http://kewaunee.com)

### Kewaunee Labway India Pvt. Ltd.

No. CA-9A, 2nd Floor, Jigani Link Road  
Jigani, Bangalore 562 106, India  
Phone: 9180-27826725 • Fax: 9180-27826724  
Web Site: [www.kewaunee.in](http://www.kewaunee.in)

### Kewaunee Labway Asia Pte. Ltd.

194 Pandan Loop, #06-22  
Pantech Business Hub, Singapore 128 383  
Phone: 65-6773-0288 • Fax: 65-6773-2322



KEWAUNEE®

Kewaunee, Alpha, BasikBench, Enterprise, Evolution, Kemresin, Kewaunee Matrix, Research Collection, Signature, and Supreme Air are registered trademarks of Kewaunee Scientific Corporation