

NAVIGATOR™

Product Overview

Top medical brands trust Cirtec Medical to accelerate time to market with vertically integrated solutions, reducing program risks. Cirtec's platform handle, Navigator, is an application-focused solution designed for

a wide range of clinical-ready catheter use cases. Engineered for simplicity and elegance, Navigator meets all necessary requirements and is ready for manufacturing.

Applications Oriented & Purpose Built

Target Applications

- Cardiovascular
- Electrophysiology
- Peripheral Vascular
- Smart Catheters
- Neurovascular

Key Product Features

- Accommodates 5Fr to 15Fr OD catheters
- Unidirectional and bidirectional steerability
- Up to 270 degrees tip deflection
- Ergonomic handle design
- Soft grip knob
- Modular design for easy customization

Flexible Configuration

Open Lumen, Flushable Catheter

- Flushable Side Port
- Proximal Seal

Closed Lumen, Smart Catheter

- Accommodates PCBA within Handle
- Configurable proximal cable and connector

Product Benefits


- Highly configurable
- Vertically integrated
- Reduced time to market
- Derisks program timelines
- Clinical ready



GUIDED WITH VERSATILITY & CONFIGURABILITY

Standards & Specifications

Up to 270 Degrees of Bi-Directional Tip Deflection



Cirtec Medical is a global leader in design, development, and manufacturing of complex Class II and III medical devices, specializing in active implant systems, interventional devices, and precision components.

Since 1987, Cirtec Medical has built extensive in-house capabilities that enable the company to serve as a single-source partner and collaborate with customers to transform their technology into innovative medical devices. With 11 global facilities and a vertically integrated approach, Cirtec Medical aims to expedite time to market, reduce risk, and accelerate revenue generation for our partners.

Learn more at CirtecMed.com