

myPatch[®]sl

Reusable, 3-channel,
Extended Wear
Holter Monitor



More
electrode
sizes



Pediatric



Neonate

Attached in under 30 seconds,
the lightweight, showerproof
myPatch[®]sl gives patients freedom
to continue daily routines without
limitations, while delivering greater
clarity and results.



Up to 14 days life on a
single battery charge

Superior in-house patient care
with no out of pocket costs.

You'll get everything you
need to start cardiac
monitoring and analysis.

- All Holter devices, electrodes and reporting software provided
- No maintenance costs or licensing fees
- Certified Cardiac Physiologist reviewed reports
- Set up, training and 24/7 customer support delivered by US based specialists

24/7 US-based customer support and direct contact with your service team provides greater certainty for your clinic and patients.

The single use, breathable fabric patch with reusable recorder allows for up to 14 days continuous cardiac monitoring in comfort, with no battery recharging patch replacement during monitoring periods.

Clinical benefits

Multi-day, multi-channel

Studies can be conducted with 2 or 3 channels, allowing for up to 14 days on 2 channels and up to 9 days on all channels.

Secure, accurate data

Greater than 99.8% average reporting accuracy, accessed and managed through BeatBox, our secure cloud-based platform.

Reusable device

Unlike single-use devices, myPatch[®]sl is reusable and only the electrodes need to be replaced.

Rapid results

With the simple to upload study data process, in as little as 24 hours you'll receive Cardiac Technician & Physiologist analyzed reports with full audit trails and multiple QA processes.

Patient benefits

Showerproof

Patients can go about daily routines without compromising data integrity or affecting their study results.

Comfort & compliance

The lightweight and compact device makes compliance easy, with many patients almost forgetting they have a device on.

3 sizes of electrodes

Adult, Pediatric and Neonate patches for patients of all ages.

Event alert function

Patients feeling symptomatic can double-tap the device to mark an event on the study,

Contact us

Cardiac Monitoring Service
1300 Quail Street, Suite 209 Newport Beach, CA 92660

T: (800)-267-5248
E: info@cmscardioscan.com

