


# Alto Holter Monitor

## What do the lights mean?

● ○ ○ **Only Green**


Green light will flash twice every 2 seconds to indicate the device is receiving a good signal and the test is underway.

After a few minutes, the green light will flash once every 2 seconds when recording starts. It will continue to flash for the test duration.




○ ● ○ **Only Yellow**

Yellow flashing light means the battery is low. The device will still record normally in this mode.




● ● ● **Green / Yellow / Red**

A series of green, yellow, and red will flash when you double-tap the device to record a symptomatic event.



○ ○ ● **Only Red**

When the red light is permanently on, please contact your doctor/clinic for support.



**Need help?**

Contact your local support team by visiting [cardiacdiagnostics.com.au](https://cardiacdiagnostics.com.au)