


Accessories

30° Angled Probe Holder and Vessel Cover



By angling the fiber optic probe at 30° into the dissolution vessel you accomplish the following:

- Reduction in hydrodynamic disturbance
- Release of air bubbles from pathlength aperture
- Reduction in ability of mirrored surface to accumulate particulates

Pump the light, not the liquid

Pion stands behind the science. Let us prove it to you. Contact us today.

Pion Inc. | 10 Cook St. | Billerica, MA 01821 | U.S.A | +1-978-528-2020 | sales@pion-inc.com | www.pion-inc.com