



VIVIT[®]
SUPPLY

PRODUCT CATALOG

vivitsupply.com



CONTENTS

- 3** Supply Chain of the Future
- 4** ASTM-Rated Medical Face Masks
- 6** Isolation and Cover Gowns
- 8** Nitrile Medical Examination Gloves
- 12** Headwear and Footwear
- 14** Sterilization Wrap

CREATING THE SUPPLY CHAIN OF THE FUTURE

An innovative factory-direct procurement model for more efficient, reliable & transparent PPE sourcing.



EXPERIENCE A NEW LEVEL OF CONTROL WITH VIVIT SUPPLY'S **FACTORY DIRECT SOURCING**

Supplier Diversification

- Multiple factory locations
- Over 20 years of PPE manufacturing experience
- Dedicated production capacity

Cost Saving

- Lowest acquisition cost
- Volume discounts
- Credit terms

Supply Chain Efficiency

- Reduce labor cost
- Streamlined logistics
- Direct factory relationship

Supply Chain Resilience

- Complete transparency
- Direct communication with factory
- Safety stock program





FACIAL PROTECTION

VIVIT Supply's facial protection line is designed to ensure clinicians are comfortable during use – both in terms of how they are protected and how they feel. Rigorously tested to the ASTM F2100 standards, VIVIT Supply's ASTM compliant face masks improve the clinician experience through industry-leading comfort and innovation. A comprehensive line for any application hospital wide, the VIVIT Supply facial protection program delivers exceptional quality and value.



**DESIGNED TO ENSURE CLINICIANS ARE
PROTECTED AND COMFORTABLE.**



PRODUCT CATALOG

Facial Protection

SKU	PRODUCT DESCRIPTION	BFE	PFE	ΔP	FLUID RESISTANCE	FLAME SPREAD	PACK FACTORS
MEDICAL FACE MASKS							
10111EM	Medical Face Mask · Ear Loops · Blue · ASTM Level 1	>95%	>80%	<5.0	80 mmHg	Class I	50 Masks/Box · 500 Masks/Case
10112EM	Medical Face Mask · Ear Loops · Yellow · ASTM Level 1	>95%	>80%	<5.0	80 mmHg	Class I	50 Masks/Box · 500 Masks/Case
10215EM	Medical Face Mask · Ear Loops · Soft · Hypoallergenic · White · ASTM Level 1	>95%	>80%	<5.0	80 mmHg	Class I	50 Masks/Box · 500 Masks/Case
11311EM	Medical Face Mask · Ear Loops · Anti-Fog Foam · Blue · ASTM Level 1	>95%	>80%	<5.0	80 mmHg	Class I	50 Masks/Box · 500 Masks/Case
20111EM	Medical Face Mask · Ear Loops · Blue · ASTM Level 2	>98%	>85%	<6.0	120 mmHg	Class I	50 Masks/Box · 500 Masks/Case
30114EM	Medical Face Mask · Ear Loops · Opal · ASTM Level 3	>98%	>85%	<6.0	160 mmHg	Class I	50 Masks/Box · 500 Masks/Case
31314EM	Medical Face Mask · Ear Loops · Anti-Fog Foam · Opal · ASTM Level 3	>98%	>85%	<6.0	160 mmHg	Class I	50 Masks/Box · 500 Masks/Case
10216PM	Medical Face Mask · Ear Loops · Pediatric · Printed · ASTM F2100 Level 1	>95%	>80%	<5.0	80 mmHg	Class I	50 Masks/Box · 500 Masks/Case
10215TM	Medical Face Mask · Ties · Soft · White · ASTM Level 1	>95%	>80%	<5.0	80 mmHg	Class I	50 Masks/Box · 300 Masks/Case
10311TM	Medical Face Mask · Ties · Blue · ASTM Level 1	>95%	>80%	<5.0	80 mmHg	Class I	50 Masks/Box · 300 Masks/Case
10317TM	Medical Face Mask · Ties · High Filtration · Grey · ASTM Level 1	>95%	>80%	<5.0	80 mmHg	Class I	50 Masks/Box · 300 Masks/Case
10333TM	Medical Face Mask · Ties · Duckbill · Green · ASTM Level 1	>95%	>80%	<5.0	80 mmHg	Class I	50 Masks/Box · 300 Masks/Case
11215TM	Medical Face Mask · Ties · Soft · Anti-fog Foam · White · ASTM Level 1	>95%	>80%	<5.0	80 mmHg	Class I	50 Masks/Box · 300 Masks/Case
11311TM	Medical Face Mask · Ties · Anti-Fog Foam · Blue · ASTM Level 1	>95%	>80%	<5.0	80 mmHg	Class I	50 Masks/Box · 300 Masks/Case
12313TM	Medical Face Mask · Ties · Anti-Fog Film · Green · ASTM Level 1	>95%	>80%	<5.0	80 mmHg	Class I	50 Masks/Box · 300 Masks/Case
13313TM	Medical Face Mask · Ties · Anti-Fog Tape · Green · ASTM Level 1	>95%	>80%	<5.0	80 mmHg	Class I	50 Masks/Box · 300 Masks/Case
20131TM	Medical Face Mask · Ties · Duckbill · Blue · ASTM Level 2	>98%	>85%	<6.0	120 mmHg	Class I	50 Masks/Box · 300 Masks/Case
20141TM	Medical Face Mask · Horizontal Ties · Blue · ASTM Level 2	>98%	>85%	<6.0	120 mmHg	Class I	50 Masks/Box · 300 Masks/Case
30215TM	Medical Face Mask · Ties · Soft · White · ASTM Level 3	>98%	>85%	<6.0	160 mmHg	Class I	50 Masks/Box · 300 Masks/Case
30314TM	Medical Face Mask · Ties · Opal · ASTM Level 3	>98%	>85%	<6.0	160 mmHg	Class I	50 Masks/Box · 300 Masks/Case
31314TM	Medical Face Mask · Ties · Anti-Fog Foam · Opal · ASTM Level 3	>98%	>85%	<6.0	160 mmHg	Class I	50 Masks/Box · 300 Masks/Case
31334TM	Medical Face Mask · Ties · Duckbill · Anti-Fog Foam · Opal · ASTM Level 3	>98%	>85%	<6.0	160 mmHg	Class I	50 Masks/Box · 300 Masks/Case
311141EMV	Medical Face Mask · Ear Loops · Anti-Glare / Anti-Fog Visor · Anti-Fog Foam · Opal · ASTM Level 3	>98%	>85%	<6.0	160 mmHg	Class I	25 Masks/Box · 100 Masks/Case
312251EMV	Medical Face Mask · Ear Loops · Anti-Glare / Anti-Fog Visor · Anti-Fog Foam · Soft · Hypoallergenic · White · ASTM Level 3	>98%	>85%	<6.0	160 mmHg	Class I	25 Masks/Box · 100 Masks/Case
313141EMV	Medical Face Mask · Ear Loops · Anti-Fog Visor · Anti-Fog Foam · Opal · ASTM Level 3	>98%	>85%	<6.0	160 mmHg	Class I	25 Masks/Box · 100 Masks/Case
201312TMV	Medical Face Mask · Ties · Duckbill · Anti-Glare / Anti-Fog Visor · Blue · ASTM Level 2	>98%	>85%	<6.0	120 mmHg	Class I	25 Masks/Box · 100 Masks/Case
313241TMV	Medical Face Mask · Ties · Anti-Fog Visor · Anti-Fog Foam · Opal · ASTM Level 3	>98%	>85%	<6.0	160 mmHg	Class I	25 Masks/Box · 100 Masks/Case
313242TMV	Medical Face Mask · Ties · Anti-Glare / Anti-Fog Visor · Anti-Fog Foam · Opal · ASTM Level 3	>98%	>85%	<6.0	160 mmHg	Class I	25 Masks/Box · 100 Masks/Case
313341TMV	Medical Face Mask · Ties · Duckbill · Anti-Glare / Anti-Fog Visor · Anti-Fog Foam · Opal · ASTM Level 3	>98%	>85%	<6.0	160 mmHg	Class I	25 Masks/Box · 100 Masks/Case
MADE IN NORTH AMERICA							
101119EM	Medical Face Mask · Ear Loops · Blue · ASTM Level 1	>95%	>80%	<5.0	80 mmHg	Class I	50 Masks/Box · 500 Masks/Case
201119EM	Medical Face Mask · Ear Loops · Blue · ASTM Level 2	>98%	>85%	<6.0	120 mmHg	Class I	50 Masks/Box · 500 Masks/Case
301149EM	Medical Face Mask · Ear Loops · Opal · ASTM Level 3	>98%	>85%	<6.0	160 mmHg	Class I	50 Masks/Box · 500 Masks/Case
313149EM	Medical Face Mask · Ear Loops · Anti-Fog Foam · Opal · ASTM Level 3	>98%	>85%	<6.0	160 mmHg	Class I	50 Masks/Box · 500 Masks/Case
102169PM	Medical Face Mask · Ear Loops · Pediatric · Printed · ASTM F2100 Level 1	>95%	>80%	<5.0	80 mmHg	Class I	50 Masks/Box · 500 Masks/Case
FACE SHIELDS & EYEWEAR							
12FS	Face Shield · Full Length · Anti-Glare / Anti-Fog Shield						10 Shields/Bag · 90 Shields/Case
22FS	Face Shield · 3/4 Length · Anti-Glare / Anti-Fog Shield						10 Shields/Bag · 90 Shields/Case
12EA	Eyewear · Assembled · Anti-Glare /Anti-Fog · Multicolored						10 Pieces/Bag · 100 Pieces/Case
12EL	Eyewear · Lens · Anti-Glare / Anti-Fog · Multicolored						100 Lenses/Box · 500 Lenses/Case
10EF	Eyewear · Frame · Reusable · Wrap Around · Black						25 Frames/Box · 200 Frames/Case





ISOLATION & COVER GOWNS

VIVIT Supply Isolation and Cover Gowns bring industry-leading innovation and protection to your gown program. Combining comfort and breathability with reliable fluid resistance, VIVIT Supply gowns are tested to global standards to provide the high level of performance clinicians expect and deserve. A wide range of neck closure systems, back styles and cuffs offer a gown suitable for any clinical use or setting.



**RELIABLE PERFORMANCE CLINICIANS
EXPECT AND DESERVE.**



PRODUCT CATALOG

Isolation Gowns

SKU	PRODUCT DESCRIPTION	PACK FACTORS
130Y	Overhead Spunbond Gown · Thumbloops with Elastic Cuff · Yellow · Universal	10 pcs/bg · 10 bg/cs · 100 pcs/cs
131Y	Overhead Spunbond Gown · Thumbloops with Elastic Cuff · Yellow · XL	10 pcs/bg · 10 bg/cs · 100 pcs/cs
220Y	SMS Cover Gown · Neck Ties · Elastic Cuff · Yellow · Universal	10 pcs/bg · 10 bg/cs · 100 pcs/cs
221Y	SMS Cover Gown · Neck Ties · Elastic Cuff · Yellow · XL	10 pcs/bg · 10 bg/cs · 100 pcs/cs
230Y	Overhead SMS Gown · Thumbloops with Elastic Cuff · Yellow · Universal	10 pcs/bg · 10 bg/cs · 100 pcs/cs
231Y	Overhead SMS Gown · Thumbloops with Elastic Cuff · Yellow · XL	10 pcs/bg · 10 bg/cs · 100 pcs/cs
240Y	Overhead AAMI 2 Gown · Open Back · Thumbloops with Elastic Cuff · Yellow · Universal	10 pcs/bg · 10 bg/cs · 100 pcs/cs
241Y	Overhead AAMI 2 Gown · Open Back · Thumbloops with Elastic Cuff · Yellow · XL	10 pcs/bg · 10 bg/cs · 100 pcs/cs
320Y	Poly-Coated Gown · Neck Ties · Knit Cuff · Yellow · Universal	10 pcs/bg · 10 bg/cs · 100 pcs/cs
321Y	Poly-Coated Gown · Neck Ties · Knit Cuff · Yellow · XL	10 pcs/bg · 10 bg/cs · 100 pcs/cs
330Y	Overhead Poly-Coated Gown · Thumbloops With Elastic Cuff · Yellow · Universal	10 pcs/bg · 10 bg/cs · 100 pcs/cs
331Y	Overhead Poly-Coated Gown · Thumbloops With Elastic Cuff · Yellow · XL	10 pcs/bg · 10 bg/cs · 100 pcs/cs
433B	Overhead AAMI 4 Gown · Open Back · Thumbloops · Blue · Universal	10pcs/bg · 10 bg/cs · 100 pcs/cs
434B	Overhead AAMI 4 Gown · Open Back · Thumbloops · Blue · XL	10pcs/bg · 10 bg/cs · 100 pcs/cs
436B	Overhead AAMI 4 Gown · Open Back · Knit Cuff · Blue · Universal	10pcs/bg · 10 bg/cs · 100 pcs/cs
437B	Overhead AAMI 4 Gown · Open Back · Knit Cuff · Blue · XL	10pcs/bg · 10 bg/cs · 100 pcs/cs
820Y	Overhead AAMI 2 Isolation Gown · Thumbloops With Elastic Cuff · Yellow · Universal	10 pcs/bg · 10 bg/cs · 100 pcs/cs
821Y	Overhead AAMI 2 Isolation Gown · Thumbloops With Elastic Cuff · Yellow · XL	10 pcs/bg · 10 bg/cs · 100 pcs/cs
824Y	Overhead AAMI 2 Isolation Gown · Thumbloops With Elastic Cuff · Yellow · Universal	10 pcs/bg · 10 bg/cs · 100 pcs/cs
825Y	Overhead AAMI 2 Isolation Gown · Thumbloops With Elastic Cuff · Yellow · XL	10 pcs/bg · 10 bg/cs · 100 pcs/cs
830Y	Overhead AAMI 3 Isolation Gown · Thumbloops With Elastic Cuff · Yellow · Universal	10 pcs/bg · 10 bg/cs · 100 pcs/cs
831Y	Overhead AAMI 3 Isolation Gown · Thumbloops With Elastic Cuff · Yellow · XL	10 pcs/bg · 10 bg/cs · 100 pcs/cs
840Y	AAMI 2 Isolation Gown · Neck Ties · Thumbloops with Elastic Cuff · Yellow · Universal	10 pcs/bg · 10 bg/cs · 100 pcs/cs
841Y	AAMI 2 Isolation Gown · Neck Ties · Thumbloops with Elastic Cuff · Yellow · XL	10 pcs/bg · 10 bg/cs · 100 pcs/cs
850Y	AAMI 3 Isolation Gown · Neck Ties · Thumbloops with Elastic Cuff · Yellow · Universal	10 pcs/bg · 10 bg/cs · 100 pcs/cs
851Y	AAMI 3 Isolation Gown · Neck Ties · Thumbloops with Elastic Cuff · Yellow · XL	10 pcs/bg · 10 bg/cs · 100 pcs/cs
860Y	AAMI 2 Isolation Gown · Tape Tab · Elastic Cuff · Yellow · Universal	10 pcs/bg · 10 bg/cs · 100 pcs/cs
861Y	AAMI 2 Isolation Gown · Tape Tab · Elastic Cuff · Yellow · XL	10 pcs/bg · 10 bg/cs · 100 pcs/cs
870Y	AAMI 3 Isolation Gown · Tape Tab · Elastic Cuff · Yellow · Universal	10 pcs/bg · 10 bg/cs · 100 pcs/cs
871Y	AAMI 3 Isolation Gown · Tape Tab · Elastic Cuff · Yellow · XL	10 pcs/bg · 10 bg/cs · 100 pcs/cs
880Y	Overhead AAMI 2 Isolation Gown · Thumbloops With Elastic Cuff · Yellow · Universal	10 pcs/bg · 10 bg/cs · 100 pcs/cs
920P	Chemotherapy Drug Tested AAMI 2 Isolation Gown · Hook & Loop Closure · Knit Cuffs · Purple · Universal	10 pcs/bg · 5 bg/cs · 50 pcs/cs
921P	Chemotherapy Drug Tested AAMI 2 Isolation Gown · Hook & Loop Closure · Knit Cuffs · Purple · XL	10 pcs/bg · 5 bg/cs · 50 pcs/cs
970B	Overhead Film Gown · Thumbloops · Blue · Universal	20 pcs/bx · 6 bx/cs · 120 pcs/cs
971B	Overhead Film Gown · Thumbloops · Blue · XL	20 pcs/bx · 6 bx/cs · 120 pcs/cs
980B	Overhead Film Gown · Thumbloops · Blue · Universal	20 pcs/bg · 6 bg/cs · 120 pcs/cs
981B	Overhead Film Gown · Thumbloops · Blue · XL	20 pcs/bg · 6 bg/cs · 120 pcs/cs





EXAMINATION GLOVES

VIVIT Supply Medical Examination Gloves use only the highest quality materials and industry leading manufacturing technology to deliver the most comfortable and protective gloves in the market. All products are tested and exceed the ASTM standards, including resistance to permeation by chemotherapy drugs and performance against viral penetration.



**DESIGNED FOR COMFORT,
PROTECTION, AND DURABILITY**



PRODUCT CATALOG

Gloves

SKU	SIZE	PRODUCT DESCRIPTION	PACK FACTORS
COMMEND			
20121	XS	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Fentanyl Rated · Indigo	XS-L gloves: 300 Gloves/Box · 10 Boxes/Case · 3000 Gloves/Case XL gloves: 250 Gloves/Box · 10 Boxes/Case · 2500 Gloves/Case
20122	S	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Fentanyl Rated · Indigo	
20123	M	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Fentanyl Rated · Indigo	
20124	L	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Fentanyl Rated · Indigo	
20125	XL	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Fentanyl Rated · Indigo	
NOVA			
20311	XS	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Fentanyl Rated · Blue	XS-L gloves: 250 Gloves/Box · 10 Boxes/Case · 2500 Gloves/Case XL gloves: 230 Gloves/Box · 10 Boxes/Case · 2300 Gloves/Case
20312	S	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Fentanyl Rated · Blue	
20313	M	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Fentanyl Rated · Blue	
20314	L	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Fentanyl Rated · Blue	
20315	XL	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Fentanyl Rated · Blue	
DEFEND			
20521	XS	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue	XS-L gloves: 200 Gloves/Box · 10 Boxes/Case · 2000 Gloves/Case XL gloves: 180 Gloves/Box · 10 Boxes/Case · 1800 Gloves/Case
20522	S	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue	
20523	M	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue	
20524	L	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue	
20525	XL	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue	
COMPLETE WITH CLEAN-DISPENSE FILM			
20541	XS	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Fentanyl Rated · Blue	250 Gloves/Box · 10 Boxes/Case · 2500 Gloves/Case
20542	S	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Fentanyl Rated · Blue	
20543	M	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Fentanyl Rated · Blue	
20544	L	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Fentanyl Rated · Blue	
20545	XL	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Fentanyl Rated · Blue	
ENHANCE			
20621	XS	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue	150 Gloves/Box · 10 Boxes/Case · 1500 Gloves/Case
20622	S	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue	
20623	M	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue	
20624	L	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue	
20625	XL	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue	
MIGHTY			
20721	XS	Nitrile Medical Examination Glove · Powder-Free · Fully Textured · Chemotherapy Drug Tested · Blue	100 Gloves/Box · 10 Boxes/Case · 1000 Gloves/Case
20722	S	Nitrile Medical Examination Glove · Powder-Free · Fully Textured · Chemotherapy Drug Tested · Blue	
20723	M	Nitrile Medical Examination Glove · Powder-Free · Fully Textured · Chemotherapy Drug Tested · Blue	
20724	L	Nitrile Medical Examination Glove · Powder-Free · Fully Textured · Chemotherapy Drug Tested · Blue	
20725	XL	Nitrile Medical Examination Glove · Powder-Free · Fully Textured · Chemotherapy Drug Tested · Blue	





PRODUCT CATALOG

Gloves

SKU	SIZE	PRODUCT DESCRIPTION	PACK FACTORS
STEALTH			
44911	XS	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Fentanyl Rated · Black	100 Gloves/Box · 10 Boxes/Case · 1000 Gloves/Case
44912	S	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Fentanyl Rated · Black	
44913	M	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Fentanyl Rated · Black	
44914	L	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Fentanyl Rated · Black	
44915	XL	Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Fentanyl Rated · Black	
LIBERTY			
40531	XS	Nitrile Medical Examination Glove · Accelerator-Free · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue	150 Gloves/Box · 10 Boxes/Case · 1500 Gloves/Case
40532	S	Nitrile Medical Examination Glove · Accelerator-Free · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue	
40533	M	Nitrile Medical Examination Glove · Accelerator-Free · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue	
40534	L	Nitrile Medical Examination Glove · Accelerator-Free · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue	
40535	XL	Nitrile Medical Examination Glove · Accelerator-Free · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue	
ESSENCE			
40011	XS	Biodegradable Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Fentanyl Rated · Violet Blue	200 Gloves/Box · 10 Boxes/Case · 2000 Gloves/Case
40012	S	Biodegradable Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Fentanyl Rated · Violet Blue	
40013	M	Biodegradable Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Fentanyl Rated · Violet Blue	
40014	L	Biodegradable Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Fentanyl Rated · Violet Blue	
40015	XL	Biodegradable Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Fentanyl Rated · Violet Blue	
BEYOND			
25911	XS	Extended Cuff Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue	100 Gloves/Box · 10 Boxes/Case · 1000 Gloves/Case
25912	S	Extended Cuff Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue	
25913	M	Extended Cuff Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue	
25914	L	Extended Cuff Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue	
25915	XL	Extended Cuff Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue	
BEYOND PRO			
26991	XS	12" Extended Cuff Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Indigo	100 Gloves/Box · 10 Boxes/Case · 1000 Gloves/Case
26992	S	12" Extended Cuff Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Indigo	
26993	M	12" Extended Cuff Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Indigo	
26994	L	12" Extended Cuff Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Indigo	
26995	XL	12" Extended Cuff Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Indigo	
EXCEL EC			
35911	XS	16" Extended Cuff Nitrile Medical Examination Glove · Powder-Free · Smooth · Chemotherapy Drug Tested · Fentanyl Rated · Blue	XS-L: 50 Gloves/Box · 10 Boxes/Case · 500 Gloves/Case XL: 45 Gloves/Box · 10 Boxes/Case · 450 Gloves/Case
35912	S	16" Extended Cuff Nitrile Medical Examination Glove · Powder-Free · Smooth · Chemotherapy Drug Tested · Fentanyl Rated · Blue	
35913	M	16" Extended Cuff Nitrile Medical Examination Glove · Powder-Free · Smooth · Chemotherapy Drug Tested · Fentanyl Rated · Blue	
35914	L	16" Extended Cuff Nitrile Medical Examination Glove · Powder-Free · Smooth · Chemotherapy Drug Tested · Fentanyl Rated · Blue	
35915	XL	16" Extended Cuff Nitrile Medical Examination Glove · Powder-Free · Smooth · Chemotherapy Drug Tested · Fentanyl Rated · Blue	





PRODUCT CATALOG

Gloves

SKU	SIZE	PRODUCT DESCRIPTION	PACK FACTORS
STERILE SINGLES			
25611	XS	Extended Cuff · Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue · Ambidextrous	100 Gloves/Box · 4 Boxes/Case · 400 Gloves/Case
25612	S	Extended Cuff · Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue · Ambidextrous	
25613	M	Extended Cuff · Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue · Ambidextrous	
25614	L	Extended Cuff · Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue · Ambidextrous	
25615	XL	Extended Cuff · Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue · Ambidextrous	
STERILE PAIRS			
25621	XS	Extended Cuff · Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue	100 Gloves/Box · 4 Boxes/Case · 400 Gloves/Case
25622	S	Extended Cuff · Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue	
25623	M	Extended Cuff · Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue	
25624	L	Extended Cuff · Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue	
25625	XL	Extended Cuff · Nitrile Medical Examination Glove · Powder-Free · Textured Fingertips · Chemotherapy Drug Tested · Blue	
STRETCH VINYL			
50711	XS	Vinyl Examination Glove · Powder-Free · Smooth · Beige	150 Gloves/Box · 10 Boxes/Case · 1500 Gloves/Case
50712	S	Vinyl Examination Glove · Powder-Free · Smooth · Beige	
50713	M	Vinyl Examination Glove · Powder-Free · Smooth · Beige	
50714	L	Vinyl Examination Glove · Powder-Free · Smooth · Beige	
50715	XL	Vinyl Examination Glove · Powder-Free · Smooth · Beige	
THRIVE			
60881	XS	Latex Examination Glove · Powder-Free · Fully Textured · Natural Color	100 Gloves/Box · 10 Boxes/Case · 1000 Gloves/Case
60882	S	Latex Examination Glove · Powder-Free · Fully Textured · Natural Color	
60883	M	Latex Examination Glove · Powder-Free · Fully Textured · Natural Color	
60884	L	Latex Examination Glove · Powder-Free · Fully Textured · Natural Color	
60885	XL	Latex Examination Glove · Powder-Free · Fully Textured · Natural Color	
THRIVE EC			
65122	S	12" Extended Cuff Latex Examination Glove · Powder-Free · Fully Textured · Chemotherapy Drug Tested · Blue	50 Gloves/Box · 10 Boxes/Case · 500 Gloves/Case
65123	M	12" Extended Cuff Latex Examination Glove · Powder-Free · Fully Textured · Chemotherapy Drug Tested · Blue	
65124	L	12" Extended Cuff Latex Examination Glove · Powder-Free · Fully Textured · Chemotherapy Drug Tested · Blue	
65125	XL	12" Extended Cuff Latex Examination Glove · Powder-Free · Fully Textured · Chemotherapy Drug Tested · Blue	





HEADWEAR & FOOTWEAR

VIVIT Supply headwear and footwear offer the ideal solution for maintaining hygienic conditions in clinical environments. Made with lightweight materials, VIVIT Supply bouffant caps are cool and breathable providing clinician comfort while preventing contamination. VIVIT Supply shoe covers come in a variety of materials offering different levels of durability and protection. Whether you're looking for basic barrier protection or protection from blood and body fluids, VIVIT Supply has a suitable option. Designed to fit all clinicians, VIVIT Supply headwear and footwear are available in a range of sizes.



**COOL & BREATHABLE WHILE
PREVENTING CONTAMINATION**



PRODUCT CATALOG

Headwear & Footwear


SKU	SIZE	PRODUCT DESCRIPTION	PACK FACTORS
1120HW	S	Bouffant Cap · SMS · Blue · Small	250 pcs/bx · 5 bx/cs · 1250 pcs/cs
1121HW	L	Bouffant Cap · SMS · Blue · Large	250 pcs/bx · 5 bx/cs · 1250 pcs/cs
1122HW	XL	Bouffant Cap · SMS · Blue · Extra-Large	100 pcs/bx · 5 bx/cs · 500 pcs/cs
1130HW	S	Bouffant Cap · SBPP · Blue · Small	250 pcs/bx · 5 bx/cs · 1250 pcs/cs
1131HW	L	Bouffant Cap · SBPP · Blue · Large	250 pcs/bx · 5 bx/cs · 1250 pcs/cs
1132HW	XL	Bouffant Cap · SBPP · Blue · Extra-Large	100 pcs/bx · 5 bx/cs · 500 pcs/cs
1150HW	S	Bouffant Cap · SBPP · White · Small	250 pcs/bx · 5 bx/cs · 1250 pcs/cs
1151HW	L	Bouffant Cap · SBPP · White · Large	250 pcs/bx · 5 bx/cs · 1250 pcs/cs
1161HW	L	Bouffant Cap · SBPP · Red · Large	250 pcs/bx · 5 bx/cs · 1250 pcs/cs
1240HW	U	Surgical Cap · Scrim · Blue · Universal	100 pcs/bx · 5 bx/cs · 500 pcs/cs
1241HW	U	Surgical Cap · SMS · Blue · Universal	100 pcs/bx · 5 bx/cs · 500 pcs/cs
1330HW	U	Surgical Hood · SMS · Tie Neck · Blue · Universal	100 pcs/bx · 3 bx/cs · 300 pcs/cs
1331HW	U	Surgical Hood · SBPP · Elastic Neck · Blue · Universal	100 pcs/bx · 3 bx/cs · 300 pcs/cs
1340HW	U	Beard Cover · Spunbond · White · Universal	100 pcs/bg · 10 bg/cs · 1000 pcs/cs
8430FW	U	Anti-Skid Shoe Cover · SMS · Blue · Universal	100 pcs/bx · 3 bx/cs · 300 pcs/cs
8431FW	XL	Anti-Skid Shoe Cover · SMS · Blue · Extra-Large	100 pcs/bx · 3 bx/cs · 300 pcs/cs
8440FW	U	Anti-Skid Shoe Cover · SBPP · Blue · Universal	100 pcs/bx · 3 bx/cs · 300 pcs/cs
8441FW	XL	Anti-Skid Shoe Cover · SBPP · Blue · Extra-Large	100 pcs/bx · 3 bx/cs · 300 pcs/cs
8450FW	U	Impervious Shoe Cover · SB/PE · Blue · Universal	100 pcs/bx · 5 bx/cs · 500 pcs/cs
8451FW	XL	Impervious Shoe Cover · SB/PE · Blue · Extra-Large	100 pcs/bx · 5 bx/cs · 500 pcs/cs
9340FW	U	Anti-Skid Boot Covers · PE Coated SMS · Knee-High · Blue · Universal	50 pcs/bx · 2 bx/cs · 100 pcs/cs
9341FW	XL	Anti-Skid Boot Covers · PE Coated SMS · Knee-High · Blue · Extra-Large	50 pcs/bx · 2 bx/cs · 100 pcs/cs





PHOENIX **STERILIZATION** **WRAP**

VIVIT Supply Phoenix Sterilization Wraps are made of high-quality SMS (spunbond-meltblown-spunbond) polypropylene, combining strength and superior filtration capabilities to provide a barrier to protect medical devices from contamination. VIVIT Supply materials offer a smooth finishing pattern to reduce friction on skin and are highly drapeable, allowing for an easy aseptic presentation.



**OFFERING BOTH STRENGTH
AND EXCEPTIONAL FILTRATION**

Available in a variety of sizes and types, there is a convenient option to meet various sterilization needs in alignment with standard hospital practices:

- VIVIT Supply Phoenix Single Layer Sterilization Wraps (S) are supplied to the customer as bulk packages of single sheets. Two sheets are then used to wrap a medical device or a collection of medical devices for sterilization.
- VIVIT Supply Phoenix Double Layer Two-color Sterilization Wraps (F2) are made of two sheets in different colours ultrasonically seamed together.

TABLE 1: DIMENSIONAL SPECIFICATIONS OF THE WRAPS

Model*	V100		V200		V300		V400		V500		V600	
Type	F2	S	F2	S	F2	S	F2	S	F2	S	F2	S
12" x 12"	•	•	•	•								
15" x 15"	•	•		•								
18" x 18"	•	•	•	•	•	•	•	•	•	•		
20" x 20"	•	•										
24" x 24"	•	•	•	•	•	•	•	•	•	•		
30" x 30"	•	•	•	•	•	•		•	•	•		
36" x 36"	•	•	•	•	•	•	•	•	•	•	•	•
40" x 40"	•	•	•	•	•	•	•	•		•	•	•
45" x 45"	•	•		•	•	•	•	•	•	•	•	•
48" x 48"	•	•		•	•	•	•	•	•	•	•	•
54" x 54"	•	•		•	•	•	•	•	•	•	•	•
60" x 60"									•			
54" x 72"	•	•	•	•	•	•	•	•	•	•	•	•
54" x 90"	•				•				•		•	

*package is color-coded by model





PRODUCT CATALOG

PHOENIX Sterilization Wrap

SKU	SIZE (IN)	WEIGHT	PRODUCT DESCRIPTION	PACK FACTORS
DOUBLE LAYER TWO COLOR • V100				
112DL	12 X 12	V100	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V100 • 12 X 12	240 pcs/bg • 2 bg/cs • 480 pcs/cs
115DL	15 X 15	V100	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V100 • 15 X 15	240 pcs/bg • 2 bg/cs • 480 pcs/cs
118DL	18 X 18	V100	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V100 • 18 X 18	240 pcs/bg • 2 bg/cs • 480 pcs/cs
120DL	20 X 20	V100	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V100 • 20 X 20	240 pcs/bg • 2 bg/cs • 480 pcs/cs
124DL	24 X 24	V100	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V100 • 24 X 24	120 pcs/bg • 2 bg/cs • 240 pcs/cs
130DL	30 X 30	V100	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V100 • 30 X 30	48 pcs/bg • 3 bg/cs • 144 pcs/cs
136DL	36 X 36	V100	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V100 • 36 X 36	48 pcs/bg • 3 bg/cs • 144 pcs/cs
140DL	40 X 40	V100	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V100 • 40 X 40	24 pcs/bg • 5 bg/cs • 120 pcs/cs
145DL	45 X 45	V100	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V100 • 45 X 45	24 pcs/bg • 5 bg/cs • 120 pcs/cs
148DL	48 X 48	V100	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V100 • 48 X 48	24 pcs/bg • 5 bg/cs • 120 pcs/cs
154DL	54 X 54	V100	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V100 • 54 X 54	24 pcs/bg • 4 bg/cs • 96 pcs/cs
190DL	54 X 90	V100	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V100 • 54 X 90	15 pcs/bg • 2 bg/cs • 30 pcs/cs
DOUBLE LAYER TWO COLOR • V200				
212DL	12 X 12	V200	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V200 • 12 X 12	240 pcs/bg • 2 bg/cs • 480 pcs/cs
218DL	18 X 18	V200	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V200 • 18 X 18	240 pcs/bg • 2 bg/cs • 480 pcs/cs
224DL	24 X 24	V200	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V200 • 24 X 24	120 pcs/bg • 2 bg/cs • 240 pcs/cs
230DL	30 X 30	V200	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V200 • 30 X 30	48 pcs/bg • 3 bg/cs • 144 pcs/cs
236DL	36 X 36	V200	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V200 • 36 X 36	48 pcs/bg • 3 bg/cs • 144 pcs/cs
240DL	40 X 40	V200	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V200 • 40 X 40	24 pcs/bg • 5 bg/cs • 120 pcs/cs
DOUBLE LAYER TWO COLOR • V300				
318DL	18 X 18	V300	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V300 • 18 X 18	240 pcs/bg • 2 bg/cs • 480 pcs/cs
324DL	24 X 24	V300	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V300 • 24 X 24	120 pcs/bg • 2 bg/cs • 240 pcs/cs
330DL	30 X 30	V300	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V300 • 30 X 30	48 pcs/bg • 3 bg/cs • 144 pcs/cs
336DL	36 X 36	V300	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V300 • 36 X 36	48 pcs/bg • 3 bg/cs • 144 pcs/cs
340DL	40 X 40	V300	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V300 • 40 X 40	24 pcs/bg • 4 bg/cs • 96 pcs/cs
345DL	45 X 45	V300	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V300 • 45 X 45	24 pcs/bg • 2 bg/cs • 48 pcs/cs
348DL	48 X 48	V300	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V300 • 48 X 48	24 pcs/bg • 4 bg/cs • 96 pcs/cs
354DL	54 X 54	V300	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V300 • 54 X 54	24 pcs/bg • 2 bg/cs • 48 pcs/cs
372DL	54 X 72	V300	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V300 • 54 X 72	24 pcs/bg • 2 bg/cs • 48 pcs/cs
390DL	54 X 90	V300	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V300 • 54 X 90	15 pcs/bg • 2 bg/cs • 30 pcs/cs
DOUBLE LAYER TWO COLOR • V400				
418DL	18 X 18	V400	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V400 • 18 X 18	144 pcs/bg • 2 bg/cs • 288 pcs/cs
424DL	24 X 24	V400	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V400 • 24 X 24	60 pcs/bg • 2 bg/cs • 120 pcs/cs
436DL	36 X 36	V400	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V400 • 36 X 36	36 pcs/bg • 2 bg/cs • 72 pcs/cs
440DL	40 X 40	V400	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V400 • 40 X 40	24 pcs/bg • 3 bg/cs • 72 pcs/cs
445DL	45 X 45	V400	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V400 • 45 X 45	24 pcs/bg • 2 bg/cs • 48 pcs/cs
448DL	48 X 48	V400	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V400 • 48 X 48	24 pcs/bg • 2 bg/cs • 48 pcs/cs
454DL	54 X 54	V400	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V400 • 54 X 54	24 pcs/bg • 2 bg/cs • 48 pcs/cs
472DL	54 X 72	V400	Sterilization Wrap • Double Layer Two Color • Purple & Blue • V400 • 54 X 72	12 pcs/bg • 2 bg/cs • 24 pcs/cs





PRODUCT CATALOG

PHOENIX Sterilization Wrap

SKU	SIZE (IN)	WEIGHT	PRODUCT DESCRIPTION	PACK FACTORS
DOUBLE LAYER TWO COLOR · V500				
518DL	18 X 18	V500	Sterilization Wrap · Double Layer Two Color · Purple & Blue · V500 · 18 X 18	144 pcs/bg · 2 bg/cs · 288 pcs/cs
524DL	24 X 24	V500	Sterilization Wrap · Double Layer Two Color · Purple & Blue · V500 · 24 X 24	60 pcs/bg · 2 bg/cs · 120 pcs/cs
530DL	30 X 30	V500	Sterilization Wrap · Double Layer Two Color · Purple & Blue · V500 · 30 X 30	60 pcs/bg · 2 bg/cs · 120 pcs/cs
536DL	36 X 36	V500	Sterilization Wrap · Double Layer Two Color · Purple & Blue · V500 · 36 X 36	36 pcs/bg · 2 bg/cs · 72 pcs/cs
545DL	45 X 45	V500	Sterilization Wrap · Double Layer Two Color · Purple & Blue · V500 · 45 X 45	24 pcs/bg · 2 bg/cs · 48 pcs/cs
548DL	48 X 48	V500	Sterilization Wrap · Double Layer Two Color · Purple & Blue · V500 · 48 X 48	24 pcs/bg · 2 bg/cs · 48 pcs/cs
554DL	54 X 54	V500	Sterilization Wrap · Double Layer Two Color · Purple & Blue · V500 · 54 X 54	24 pcs/bg · 2 bg/cs · 48 pcs/cs
560DL	60 X 60	V500	Sterilization Wrap · Double Layer Two Color · Purple & Blue · V500 · 60 X 60	12 pcs/bg · 2 bg/cs · 24 pcs/cs
572DL	54 X 72	V500	Sterilization Wrap · Double Layer Two Color · Purple & Blue · V500 · 54 X 72	12 pcs/bg · 2 bg/cs · 24 pcs/cs
590DL	54 X 90	V500	Sterilization Wrap · Double Layer Two Color · Purple & Blue · V500 · 54 X 90	12 pcs/bg · 2 bg/cs · 24 pcs/cs
DOUBLE LAYER TWO COLOR · V600				
636DL	36 X 36	V600	Sterilization Wrap · Double Layer Two Color · Purple & Blue · V600 · 36 X 36	36 pcs/bg · 2 bg/cs · 72 pcs/cs
640DL	40 X 40	V600	Sterilization Wrap · Double Layer Two Color · Purple & Blue · V600 · 40 X 40	24 pcs/bg · 2 bg/cs · 48 pcs/cs
645DL	45 X 45	V600	Sterilization Wrap · Double Layer Two Color · Purple & Blue · V600 · 45 X 45	24 pcs/bg · 2 bg/cs · 48 pcs/cs
648DL	48 X 48	V600	Sterilization Wrap · Double Layer Two Color · Purple & Blue · V600 · 48 X 48	24 pcs/bg · 2 bg/cs · 48 pcs/cs
654DL	54 X 54	V600	Sterilization Wrap · Double Layer Two Color · Purple & Blue · V600 · 54 X 54	24 pcs/bg · 2 bg/cs · 48 pcs/cs
672DL	54 X 72	V600	Sterilization Wrap · Double Layer Two Color · Purple & Blue · V600 · 54 X 72	24 pcs/bg · 1 bg/cs · 24 pcs/cs
690DL	54 X 90	V600	Sterilization Wrap · Double Layer Two Color · Purple & Blue · V600 · 54 X 90	12 pcs/bg · 2 bg/cs · 24 pcs/cs
SINGLE LAYER · V100				
112SL	12 X 12	V100	Sterilization Wrap · Single Layer · Purple · V100 · 12 X 12	100 pcs/bg · 10 bg/cs · 1000 pcs/cs
115SL	15 X 15	V100	Sterilization Wrap · Single Layer · Purple · V100 · 15 X 15	100 pcs/bg · 10 bg/cs · 1000 pcs/cs
118SL	18 X 18	V100	Sterilization Wrap · Single Layer · Purple · V100 · 18 X 18	100 pcs/bg · 10 bg/cs · 1000 pcs/cs
120SL	20 X 20	V100	Sterilization Wrap · Single Layer · Purple · V100 · 20 X 20	100 pcs/bg · 10 bg/cs · 1000 pcs/cs
124SL	24 X 24	V100	Sterilization Wrap · Single Layer · Purple · V100 · 24 X 24	100 pcs/bg · 5 bg/cs · 500 pcs/cs
130SL	30 X 30	V100	Sterilization Wrap · Single Layer · Purple · V100 · 30 X 30	150 pcs/bg · 2 bg/cs · 300 pcs/cs
136SL	36 X 36	V100	Sterilization Wrap · Single Layer · Purple · V100 · 36 X 36	150 pcs/bg · 2 bg/cs · 300 pcs/cs
140SL	40 X 40	V100	Sterilization Wrap · Single Layer · Purple · V100 · 40 X 40	50 pcs/bg · 5 bg/cs · 250 pcs/cs
145SL	45 X 45	V100	Sterilization Wrap · Single Layer · Purple · V100 · 45 X 45	50 pcs/bg · 5 bg/cs · 250 pcs/cs
148SL	48 X 48	V100	Sterilization Wrap · Single Layer · Purple · V100 · 48 X 48	50 pcs/bg · 5 bg/cs · 250 pcs/cs
154SL	54 X 54	V100	Sterilization Wrap · Single Layer · Purple · V100 · 54 X 54	50 pcs/bg · 2 bg/cs · 100 pcs/cs
172SL	54 X 72	V100	Sterilization Wrap · Single Layer · Purple · V100 · 54 X 72	50 pcs/bg · 2 bg/cs · 100 pcs/cs
SINGLE LAYER · V200				
212SL	12 X 12	V200	Sterilization Wrap · Single Layer · Purple · V200 · 12 X 12	500 pcs/bg · 2 bg/cs · 1000 pcs/cs
215SL	15 X 15	V200	Sterilization Wrap · Single Layer · Purple · V200 · 15 X 15	500 pcs/bg · 2 bg/cs · 1000 pcs/cs
218SL	18 X 18	V200	Sterilization Wrap · Single Layer · Purple · V200 · 18 X 18	500 pcs/bg · 2 bg/cs · 1000 pcs/cs
224SL	24 X 24	V200	Sterilization Wrap · Single Layer · Purple · V200 · 24 X 24	250 pcs/bg · 2 bg/cs · 500 pcs/cs
230SL	30 X 30	V200	Sterilization Wrap · Single Layer · Purple · V200 · 30 X 30	100 pcs/bg · 3 bg/cs · 300 pcs/cs
236SL	36 X 36	V200	Sterilization Wrap · Single Layer · Purple · V200 · 36 X 36	100 pcs/bg · 3 bg/cs · 300 pcs/cs
240SL	40 X 40	V200	Sterilization Wrap · Single Layer · Purple · V200 · 40 X 40	50 pcs/bg · 5 bg/cs · 250 pcs/cs
248SL	48 X 48	V200	Sterilization Wrap · Single Layer · Purple · V200 · 48 X 48	50 pcs/bg · 2 bg/cs · 100 pcs/cs
272SL	54 X 72	V200	Sterilization Wrap · Single Layer · Purple · V200 · 54 X 72	50 pcs/bg · 2 bg/cs · 100 pcs/cs





PRODUCT CATALOG

PHOENIX Sterilization Wrap

SKU	SIZE (IN)	WEIGHT	PRODUCT DESCRIPTION	PACK FACTORS
SINGLE LAYER • V300				
318SL	18 X 18	V300	Sterilization Wrap • Single Layer • Purple • V300 • 18 X 18	500 pcs/bg • 2 bg/cs • 1000 pcs/cs
324SL	24 X 24	V300	Sterilization Wrap • Single Layer • Purple • V300 • 24 X 24	250 pcs/bg • 2 bg/cs • 500 pcs/cs
330SL	30 X 30	V300	Sterilization Wrap • Single Layer • Purple • V300 • 30 X 30	100 pcs/bg • 3 bg/cs • 300 pcs/cs
336SL	36 X 36	V300	Sterilization Wrap • Single Layer • Purple • V300 • 36 X 36	100 pcs/bg • 3 bg/cs • 300 pcs/cs
340SL	40 X 40	V300	Sterilization Wrap • Single Layer • Purple • V300 • 40 X 40	50 pcs/bg • 5 bg/cs • 250 pcs/cs
345SL	45 X 45	V300	Sterilization Wrap • Single Layer • Purple • V300 • 45 X 45	50 pcs/bg • 2 bg/cs • 100 pcs/cs
348SL	48 X 48	V300	Sterilization Wrap • Single Layer • Purple • V300 • 48 X 48	50 pcs/bg • 2 bg/cs • 100 pcs/cs
354SL	54 X 54	V300	Sterilization Wrap • Single Layer • Purple • V300 • 54 X 54	50 pcs/bg • 2 bg/cs • 100 pcs/cs
372SL	54 X 72	V300	Sterilization Wrap • Single Layer • Purple • V300 • 54 X 72	50 pcs/bg • 2 bg/cs • 100 pcs/cs
SINGLE LAYER • V400				
418SL	18 X 18	V400	Sterilization Wrap • Single Layer • Purple • V400 • 18 X 18	300 pcs/bg • 2 bg/cs • 600 pcs/cs
424SL	24 X 24	V400	Sterilization Wrap • Single Layer • Purple • V400 • 24 X 24	125 pcs/bg • 2 bg/cs • 250 pcs/cs
430SL	30 X 30	V400	Sterilization Wrap • Single Layer • Purple • V400 • 30 X 30	50 pcs/bg • 4 bg/cs • 200 pcs/cs
436SL	36 X 36	V400	Sterilization Wrap • Single Layer • Purple • V400 • 36 X 36	150 pcs/bg • 2 bg/cs • 300 pcs/cs
440SL	40 X 40	V400	Sterilization Wrap • Single Layer • Purple • V400 • 40 X 40	50 pcs/bg • 3 bg/cs • 150 pcs/cs
445SL	45 X 45	V400	Sterilization Wrap • Single Layer • Purple • V400 • 45 X 45	50 pcs/bg • 2 bg/cs • 100 pcs/cs
448SL	48 X 48	V400	Sterilization Wrap • Single Layer • Purple • V400 • 48 X 48	50 pcs/bg • 1 bg/cs • 50 pcs/cs
454SL	54 X 54	V400	Sterilization Wrap • Single Layer • Purple • V400 • 54 X 54	50 pcs/bg • 1 bg/cs • 50 pcs/cs
SINGLE LAYER • V500				
518SL	18 X 18	V500	Sterilization Wrap • Single Layer • Purple • V500 • 18 X 18	300 pcs/bg • 2 bg/cs • 600 pcs/cs
524SL	24 X 24	V500	Sterilization Wrap • Single Layer • Purple • V500 • 24 X 24	125 pcs/bg • 2 bg/cs • 250 pcs/cs
530SL	30 X 30	V500	Sterilization Wrap • Single Layer • Purple • V500 • 30 X 30	125 pcs/bg • 2 bg/cs • 250 pcs/cs
536SL	36 X 36	V500	Sterilization Wrap • Single Layer • Purple • V500 • 36 X 36	75 pcs/bg • 2 bg/cs • 150 pcs/cs
545SL	45 X 45	V500	Sterilization Wrap • Single Layer • Purple • V500 • 45 X 45	50 pcs/bg • 2 bg/cs • 100 pcs/cs
548SL	48 X 48	V500	Sterilization Wrap • Single Layer • Purple • V500 • 48 X 48	50 pcs/bg • 1 bg/cs • 50 pcs/cs
554SL	54 X 54	V500	Sterilization Wrap • Single Layer • Purple • V500 • 54 X 54	50 pcs/bg • 1 bg/cs • 50 pcs/cs
572SL	54 X 72	V500	Sterilization Wrap • Single Layer • Purple • V500 • 54 X 72	50 pcs/bg • 1 bg/cs • 50 pcs/cs
SINGLE LAYER • V600				
636SL	36 X 36	V600	Sterilization Wrap • Single Layer • Purple • V500 • 36 X 36	72 pcs/bg • 2 bg/cs • 144 pcs/cs
640SL	40 X 40	V600	Sterilization Wrap • Single Layer • Purple • V500 • 40 X 40	48 pcs/bg • 2 bg/cs • 96 pcs/cs
645SL	45 X 45	V600	Sterilization Wrap • Single Layer • Purple • V500 • 45 X 45	48 pcs/bg • 2 bg/cs • 96 pcs/cs
648SL	48 X 48	V600	Sterilization Wrap • Single Layer • Purple • V500 • 48 X 48	48 pcs/bg • 1 bg/cs • 48 pcs/cs
654SL	54 X 54	V600	Sterilization Wrap • Single Layer • Purple • V500 • 54 X 54	48 pcs/bg • 1 bg/cs • 48 pcs/cs
672SL	54 X 72	V600	Sterilization Wrap • Single Layer • Purple • V500 • 54 X 72	24 pcs/bg • 2 bg/cs • 48 pcs/cs
ACCESSORIES				
TLW25L20	20 X 25		Tray Liner • White • 20 X 25	45 pcs/bg • 10 bg/cs • 450 pcs/cs





VIVIT Supply LLC.
200 Continental Dr.
Newark, DE 19713 USA

P +1.302.485.0580
W vivitsupply.com
E info@vivitsupply.com