

NITRILE MEDICAL EXAMINATION GLOVES

**POWDER-FREE • TEXTURED FINGERTIPS • CHEMOTHERAPY DRUG TESTED
FENTANYL RATED • INDIGO**



2012X

Experience the perfect balance of comfort and protection with the VIVIT Supply Commend Nitrile Examination Gloves. These gloves are ideal for precision tasks that require tactile sensitivity. Utilizing cutting-edge nitrile technology, they offer a lightweight design with a snug and comfortable fit. Stand out with this advanced option that provides a high level of protection for various healthcare settings.

PROTECTION STANDARDS

ASTM D6319 Standard Specification for Nitrile Examination Gloves for Medical Application.

ASTM D6978 Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs.

ASTM F1671 Standard Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Blood-Borne Pathogens Using Phi-X174 Bacteriophage Penetration as a Test System.

SKU	SIZE	GLOVES/BOX	BOXES/CASE	GLOVES/CASE
20121	XS	300	10	3000
20122	S	300	10	3000
20123	M	300	10	3000
20124	L	300	10	3000
20125	XL	250	10	2500

① Low-Profile Nitrile Technology

- Feel confident in the ability to maintain tactile sensitivity for intricate and precise tasks while not compromising material strength.

② Motion Enhancing Fit

- Natural fit and feel conform to the user's hands to provide exceptional dexterity.

③ Full Micro Texturing and Elevated Fingertip Texturing

- Enhances grip in wet and dry conditions when handling small instruments and sensitive tasks.

