

NITRILE MEDICAL EXAMINATION GLOVES

**POWDER-FREE • TEXTURED FINGERTIPS • CHEMOTHERAPY DRUG TESTED
FENTANYL RATED • BLUE**



2031X

Discover the balance between strength, accuracy, and comfort through the utilization of lightweight nitrile technology. The VIVIT Supply Nova glove delivers the wearer with the tactile sensitivity essential for executing intricate tasks, alongside the superior level of protection necessary for healthcare providers.

PROTECTION STANDARDS

ASTM D6319 Standard Specification for Nitrile Examination Gloves for Medical Application.

ASTM D6978 Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs.

ASTM F1671 Standard Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Blood-Borne Pathogens Using Phi-X174 Bacteriophage Penetration as a Test System.

| SKU | SIZE | GLOVES/BOX | BOXES/CASE | GLOVES/CASE |
|-------|------|------------|------------|-------------|
| 20311 | XS | 250 | 10 | 2500 |
| 20312 | S | 250 | 10 | 2500 |
| 20313 | M | 250 | 10 | 2500 |
| 20314 | L | 250 | 10 | 2500 |
| 20315 | XL | 230 | 10 | 2300 |

① Low-Profile Nitrile Technology

- Feel confident in the ability to maintain tactile sensitivity for intricate and precise tasks while not compromising material strength.

② Motion Enhancing Fit

- Natural fit and feel conform to the user's hands to provide exceptional dexterity.

③ Full Micro Texturing and Elevated Fingertip Texturing

- Enhances grip in wet and dry conditions when handling small instruments and sensitive tasks.

