

NITRILE MEDICAL EXAMINATION GLOVES

POWDER-FREE • TEXTURED FINGERTIPS • CHEMOTHERAPY DRUG TESTED • BLUE



2054X

Experience maximum comfort, excellent dexterity, and reliable protection with the VIVIT Supply Complete Nitrile Examination Gloves. Trusted by clinicians, this glove can be used throughout a variety of healthcare settings and procedures. Using the latest nitrile technology, this glove offers a flexible yet durable feel while exceeding medical grade standards for hand protection and safety requirements.

PROTECTION STANDARDS

ASTM D6319 Standard Specification for Nitrile Examination Gloves for Medical Application.

ASTM D6978 Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs.

ASTM F1671 Standard Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Blood-Borne Pathogens Using Phi-X174 Bacteriophage Penetration as a Test System.

SKU	SIZE	GLOVES/BOX	BOXES/CASE	GLOVES/CASE
20541	XS	250	10	2500
20542	S	250	10	2500
20543	M	250	10	2500
20544	L	250	10	2500
20545	XL	250	10	2500

① Low-Profile Nitrile Technology

- Feel confident in the ability to maintain tactile sensitivity for intricate and precise tasks while not compromising material strength.

② Motion Enhancing Fit

- Natural fit and feel conform to the user's hands to provide exceptional dexterity.

③ Full Micro Texturing and Elevated Fingertip Texturing

- Enhances grip in wet and dry conditions when handling small instruments and sensitive tasks.

④ Clean Dispense Packaging

- The clean dispense film covers the box opening to control glove dispensing and keep debris out.

