

NITRILE MEDICAL EXAMINATION GLOVES

POWDER-FREE • TEXTURED FINGERTIPS • CHEMOTHERAPY DRUG TESTED • BLUE



2062X

When elevated protection is required, and comfort and proper fit is of utmost importance, the VIVIT Supply Enhance Nitrile Examination Gloves are an ideal choice for various medical applications. Its versatile design provides a combination of comfort, protection, and durability, making this an industry-leading glove for healthcare professionals.

PROTECTION STANDARDS

ASTM D6319 Standard Specification for Nitrile Examination Gloves for Medical Application.

ASTM D6978 Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs.

ASTM F1671 Standard Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Blood-Borne Pathogens Using Phi-X174 Bacteriophage Penetration as a Test System.

SKU	SIZE	GLOVES/BOX	BOXES/CASE	GLOVES/CASE
20621	XS	150	10	1500
20622	S	150	10	1500
20623	M	150	10	1500
20624	L	150	10	1500
20625	XL	150	10	1500
20626	XXL	130	10	1300

① Low-Profile Nitrile Technology

- Feel confident in the ability to maintain tactile sensitivity for intricate and precise tasks while not compromising material strength.

② Motion Enhancing Fit

- Natural fit and feel conform to the user's hands to provide exceptional dexterity.

③ Full Micro Texturing and Elevated Fingertip Texturing

- Enhances grip in wet and dry conditions when handling small instruments and sensitive tasks.

