

## NITRILE MEDICAL EXAMINATION GLOVES

POWDER-FREE • FULLY TEXTURED • CHEMOTHERAPY DRUG TESTED • BLUE



**2072X**

Ensure reliable defense in tough situations with the VIVIT Supply Mighty Nitrile Examination Glove. Crafted from high quality nitrile material, exceptional strength and durability are prioritized for dependable protection during procedures. Enhanced comfort, flexibility, and tactile sensitivity make this an essential glove for healthcare professionals.

### PROTECTION STANDARDS

**ASTM D6319** Standard Specification for Nitrile Examination Gloves for Medical Application.

**ASTM D6978** Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs.

**ASTM F1671** Standard Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Blood-Borne Pathogens Using Phi-X174 Bacteriophage Penetration as a Test System.

SKU	SIZE	GLOVES/BOX	BOXES/CASE	GLOVES/CASE
20721	XS	100	10	1000
20722	S	100	10	1000
20723	M	100	10	1000
20724	L	100	10	1000
20725	XL	100	10	1000

### ① Durable Protection

- Vitalizes high strength yet responsive nitrile material ideal for increased protection in higher-risk tasks that require a thicker glove.

### ② Motion Enhancing Fit

- Natural fit and feel conform to the user's hands to provide exceptional dexterity.

### ③ Fully Textured

- Enhances grip in wet and dry conditions when handling small instruments and sensitive tasks.

