

EXTENDED CUFF NITRILE MEDICAL EXAMINATION GLOVES

TEXTURED FINGERTIPS • CHEMOTHERAPY DRUG TESTED • POWDER FREE • BLUE



2562X (PAIRS)

VIVIT Supply Sterile Nitrile Examination Gloves ensure a high standard of protection and comfort, ideal for various medical or laboratory applications. Created using a premium nitrile formulation, this glove is sterilized and packaged for use during procedures where maintaining a germ-free environment is crucial. For reliable barrier protection against contaminants, this formula achieves both safety and ease of use with its durable yet flexible fit.

PROTECTION STANDARDS

ASTM D6319 Standard Specification for Nitrile Examination Gloves for Medical Application.

ASTM D6978 Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs.

ASTM F1671 Standard Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Blood-Borne Pathogens Using Phi-X174 Bacteriophage Penetration as a Test System.

SKU	SIZE	PAIRS/BOX	BOXES/CASE	GLOVES/CASE
25621	XS	50	4	400
25622	S	50	4	400
25623	M	50	4	400
25624	L	50	4	400
25625	XL	50	4	400

① Durable Protection

- Vitalizes high strength yet responsive nitrile material ideal for increased protection in higher-risk tasks that require a thicker glove.

② Motion Enhancing Fit

- Natural fit and feel conform to the user's hands to provide exceptional dexterity.

③ Full Micro Texturing and Elevated Fingertip Texturing

- Enhances grip in wet and dry conditions when handling small instruments and sensitive tasks.

④ Extended Cuff

- Longer 12" cuff provides confidence in protection through enhanced wrist coverage.

