

# BIODEGRADABLE NITRILE MEDICAL EXAMINATION GLOVES

**POWDER-FREE • TEXTURED FINGERTIPS • CHEMOTHERAPY DRUG TESTED  
FENTANYL RATED • VIOLET BLUE**


**4001X**

Step into a more sustainable future with the VIVIT Supply Essence Biodegradable Nitrile Examination Gloves. Formulated to biodegrade 81% over 490 days\*, this glove still maintains the same reliable protection and comfort that clinicians need. With the VIVIT Supply Essence there is no need to compromise quality for the health of our planet.

## PROTECTION STANDARDS

**ASTM D6319** Standard Specification for Nitrile Examination Gloves for Medical Application.

**ASTM D6978** Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs.

**ASTM F1671** Standard Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Blood-Borne Pathogens Using Phi-X174 Bacteriophage Penetration as a Test System.

SKU	SIZE	GLOVES/BOX	BOXES/CASE	GLOVES/CASE
40011	XS	200	10	2000
40012	S	200	10	2000
40013	M	200	10	2000
40014	L	200	10	2000
40015	XL	200	10	2000

### ① Low-Profile Nitrile Technology

- Feel confident in the ability to maintain tactile sensitivity for intricate and precise tasks while not compromising material strength.

### ② Motion Enhancing Fit

- Natural fit and feel conform to the user's hands to provide exceptional dexterity.

### ③ Full Micro Texturing and Elevated Fingertip Texturing

- Enhances grip in wet and dry conditions when handling small instruments and sensitive tasks.

### ④ Innovative Biodegradable Formula

- Designed to attract microbes present in landfills to promote enzymatic biodegradation, reducing waste accumulation.

\*Tested to ASTM D5526 and D5511 for biodegradability.

