

NITRILE MEDICAL EXAMINATION GLOVES

POWDER-FREE • TEXTURED FINGERTIPS • CHEMOTHERAPY DRUG TESTED • BLUE



4053X

Reduce the potential for Type IV dermatitis reactions without compromising comfort and fit using the VIVIT Supply Liberty Nitrile Examination Gloves. Uniquely designed without the use of chemical accelerators known to cause irritation, this glove offers an additional layer of protection against material sensitivities while maintaining industry-leading protection against contaminants and harmful substances.

PROTECTION STANDARDS

ASTM D6319 Standard Specification for Nitrile Examination Gloves for Medical Application.

ASTM D6978 Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs.

ASTM F1671 Standard Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Blood-Borne Pathogens Using Phi-X174 Bacteriophage Penetration as a Test System.

SKU	SIZE	GLOVES/BOX	BOXES/CASE	GLOVES/CASE
40531	XS	150	10	1500
40532	S	150	10	1500
40533	M	150	10	1500
40534	L	150	10	1500
40535	XL	150	10	1500

① Low-Profile Nitrile Technology

- Feel confident in the ability to maintain tactile sensitivity for intricate and precise tasks while not compromising material strength.

② Motion Enhancing Fit

- Natural fit and feel conform to the user's hands to provide exceptional dexterity.

③ Full Micro Texturing and Elevated Fingertip Texturing

- Enhances grip in wet and dry conditions when handling small instruments and sensitive tasks.

④ Dermatitis Prevention

- Reduced risk of Type IV skin irritations through the removal of accelerators in the manufacturing process.

* * * * *