

EXTENDED CUFF LATEX EXAMINATION GLOVES

POWDER-FREE • FULLY TEXTURED • CHEMOTHERAPY DRUG TESTED • BLUE



6512X

Adapting to the contours of your hand for enhanced comfort and tactile sensitivity, the VIVIT Supply Thrive EC Latex Examination Glove offers a reliable and hygienic barrier across a range of clinical applications. Offering an increased area of protection and durable construction, this extended cuff glove provides a convenient solution for short-term tasks requiring hand protection in both wet and dry conditions.

PROTECTION STANDARDS

ASTM D3578 Standard specification for rubber examination gloves.

ASTM D6978 Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs.

ASTM F1671 Standard Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Blood-Borne Pathogens Using Phi-X174 Bacteriophage Penetration as a Test System.

SKU	SIZE	GLOVES/BOX	BOXES/CASE	GLOVES/CASE
65122	S	50	10	500
65123	M	50	10	500
65124	L	50	10	500
65125	XL	50	10	500

① Polymer-coated Latex

- Provides various benefits including seamless donning and doffing.

② Motion Enhancing Fit

- Natural fit and feel conform to the user's hand to provide exceptional dexterity.

③ Fully Textured

- Enhances grip in wet and dry conditions when handling small instruments and sensitive tasks.

④ Extended Cuff

- Longer 12" cuff provides confidence in protection through enhanced wrist coverage.

