

MEDICAL FACE MASK

EARLOOPS • YELLOW

ASTM LEVEL 1



10112EM

For use in isolation cases or situations with a low risk of fluid exposure, clinicians can rely on this ASTM Level 1 mask to provide comfort and breathability while acting as a barrier against contaminants. Grab this mask for patient procedures including dressing changes or line changes.

BFE	PFE	DELTA P $\Delta P(\text{mmH}_2\text{O}/\text{cm}^2)$	FLUID RESISTANCE	FLAME SPREAD
> 95%	> 80%	< 5.0	80 mm Hg	Class 1
MASKS/BOX		BOXES/CASE		MASKS/CASE
50		10		500

① Ultra-soft Earloop

- Earloop material is soft and smooth for maximum comfort.
- Round design makes this mask one of the most comfortable options available.

② High Quality Nosepiece

- Aluminum nosepiece forms smoothly across the nose and cheekbones.
- Remains securely in place and reduces fogging of eyewear.

③ Responsive Inner Layer Material

- Padded feel for ultimate comfort.
- Flexible material that moves with you.



VIVIT Supply LLC

200 Continental Dr., Suite 401
Newark, DE 19713 USA

W vivitsupply.com

V004