



MEDICAL FACE MASK

TIES • SOFT • WHITE

ASTM LEVEL 1



10215TM

Acting as a physical barrier against contaminants, this ASTM Level 1 mask offers both fluid resistance and filtration efficiency to keep clinicians protected during surgical procedures. The premium inner layer — a soft, irritant-free material — makes this the perfect mask for individuals with skin sensitivities.

BFE	PFE	DELTA P $\Delta P(\text{mmH}_2\text{O}/\text{cm}^2)$	FLUID RESISTANCE	FLAME SPREAD
> 95%	> 80%	< 5.0	80 mm Hg	Class 1
MASKS/BOX		BOXES/CASE		MASKS/CASE
50		6		300

① Premium Tie Material and Design

- Tie material that is comfortable on the skin throughout extended wear.
- Long ties create a flexible fit while tying over top of a bouffant or surgical cap.

② High Quality Nosepiece

- Aluminum nosepiece forms smoothly across the nose and cheekbones.
- Remains securely in place and reduces fogging of eyewear.

③ Soft Irritant Free Inner Layer

- Inner layer has a silky finish that is delicate yet comfortable, making it ideal for those with skin sensitivities.

