



# MEDICAL FACE MASK

## EARLOOPS • PEDIATRIC • PRINTED

### ASTM F2100 LEVEL 1



10216PM

Designed to create an effective barrier against germs, fluids, and other pathogens, VIVIT Supply's ASTM Level 1 Pediatric Medical Face Masks offer an optimal fit for children to maximize protection. Crafted using vibrant and soft materials, comfort and breathability are prioritized to encourage consistent wear. As an ideal choice for limiting the spread of respiratory symptoms, grab this mask for reliable safeguarding of children's health.

BFE	PFE	DELTA P ΔP(mmH <sub>2</sub> O/cm <sup>2</sup> )	FLUID RESISTANCE	FLAME SPREAD
> 95%	> 80%	< 5.0	80 mm Hg	Class 1
MASKS/BOX		BOXES/CASE		MASKS/CASE
50		10		500

① **Ultra-soft Earloop**

- Earloop material is soft and smooth for maximum comfort.
- Round design makes this mask one of the most comfortable options available.

② **High Quality Nosepiece**

- Aluminum nosepiece forms smoothly across the nose and cheekbones.
- Remains securely in place and reduces fogging of eyewear.

③ **Soft Irritant-Free Inner Layer Material**

- Inner layer has a silky finish that is delicate yet comfortable, making it ideal for those with skin sensitivities.

