



# MEDICAL FACE MASK TIES • BLUE ASTM LEVEL 1



**10311TM**

For use in environments with minimal fluid exposure, this ASTM Level 1 mask establishes a physical barrier between clinicians and contaminants from their surrounding environment. Highly breathable materials create a lightweight feel keeping clinicians comfortable throughout surgical procedures.

BFE	PFE	DELTA P $\Delta P(\text{mmH}_2\text{O}/\text{cm}^2)$	FLUID RESISTANCE	FLAME SPREAD
> 95%	> 80%	< 5.0	80 mm Hg	Class 1
MASKS/BOX		BOXES/CASE		MASKS/CASE
50		6		300

## ① Premium Tie Material and Design

- Tie material that is comfortable on the skin throughout extended wear.
- Long ties create a flexible fit while tying over top of a bouffant or surgical cap.

## ② High Quality Nosepiece

- Aluminum nosepiece forms smoothly across the nose and cheekbones.
- Remains securely in place and reduces fogging of eyewear.

## ③ Smooth Naturally Breathable Inner Layer

- Provides a light and cool feel to the mask, reducing chances of irritation throughout extended wear.
- Inner layer maintains its structure, allowing the mouth to feel open and airy.



Not made with  
natural rubber latex



\* \* \* \* \*

**VIVIT Supply LLC**

200 Continental Dr., Suite 401  
Newark, DE 19713 USA

**W** [vivitsupply.com](http://vivitsupply.com)

V004