



# MEDICAL FACE MASK

## TIES • SOFT • ANTI-FOG FOAM • WHITE

### ASTM LEVEL 1



**11215TM**

Clinicians can rely on this ASTM Level 1 mask to provide comfort and breathability while acting as a barrier against contaminants. The addition of a soft anti-fog foam strip reduces the fogging of eyewear, ensuring a clear line of sight during surgical procedures. The premium inner layer—a soft, irritant-free material—also makes this the perfect mask for individuals with skin sensitivities.

BFE	PFE	DELTA P $\Delta P(\text{mmH}_2\text{O}/\text{cm}^2)$	FLUID RESISTANCE	FLAME SPREAD
> 95%	> 80%	< 5.0	80 mm Hg	Class 1
MASKS/BOX		BOXES/CASE		MASKS/CASE
50		6		300

#### ① Premium Tie Material and Design

- Tie material that is comfortable on the skin throughout extended wear.
- Long ties create a flexible fit while tying over top of a bouffant or surgical cap.

#### ② High Quality Nosepiece with Ultra Soft Foam

- Aluminum nosepiece and soft foam strip work together to form smoothly across the nose and cheekbones.
- Remains securely in place to create a custom and secure seal providing enhanced resistance to fogging of eyewear and a cushion-like feel to the mask.

#### ③ Soft Irritant Free Inner Layer

- Inner layer has a silky finish that is delicate yet comfortable, making it ideal for those with skin sensitivities.



Not made with  
natural rubber latex



\* \* \* \* \*

**VIVIT Supply LLC**

200 Continental Dr., Suite 401  
Newark, DE 19713 USA

**W** [vivitsupply.com](http://vivitsupply.com)

V004