

MEDICAL FACE MASK

EARLOOPS • ANTI-FOG FOAM • BLUE

ASTM LEVEL 1



11311EM

Clinicians can rely on this ASTM Level 1 mask to provide comfort and breathability while acting as a barrier against contaminants. The addition of a soft anti-fog foam strip reduces the fogging of eyewear, ensuring a clear line of sight during surgical procedures.

BFE	PFE	DELTA P $\Delta P(\text{mmH}_2\text{O}/\text{cm}^2)$	FLUID RESISTANCE	FLAME SPREAD
> 95%	> 80%	< 5.0	80 mm Hg	Class 1

MASKS/BOX	BOXES/CASE	MASKS/CASE
50	10	500

① Ultra Soft Earloop

- Earloop material is soft and smooth for maximum comfort.
- Round design makes this mask one of the most comfortable options available.

② High Quality Nosepiece with Ultra Soft Foam

- Aluminum nosepiece and soft foam strip work together to form smoothly across the nose and cheekbones.
- Remains securely in place to create a custom and secure seal providing enhanced resistance to fogging of eyewear and a cushion-like feel to the mask.

③ Responsive Inner Layer Material

- Padded feel for ultimate comfort.
- Flexible material that moves with you.



Not made with
natural rubber latex



* * * * *

VIVIT Supply LLC

200 Continental Dr., Suite 401
Newark, DE 19713 USA

W vivitsupply.com

V002