



MEDICAL FACE MASK

TIES • ANTI-FOG FILM • GREEN

ASTM LEVEL 1



12313TM

This ASTM Level 1 mask offers standard fluid resistance and filtration efficiency to keep clinicians protected during surgical procedures with a low risk of fluid exposure. Soft plastic film lines the inside and outside of the mask and is designed to cling to the skin to form a seal and act as a vapor barrier against moisture to prevent the fogging of eyewear.

BFE	PFE	DELTA P $\Delta P(\text{mmH}_2\text{O}/\text{cm}^2)$	FLUID RESISTANCE	FLAME SPREAD
> 95%	> 80%	< 5.0	80 mm Hg	Class 1
MASKS/BOX		BOXES/CASE		MASKS/CASE
50		6		300

① Premium Tie Material and Design

- Tie material that is comfortable on the skin throughout extended wear.
- Long ties create a flexible fit while tying over top of a bouffant or surgical cap.

② High Quality Nosepiece with Anti-Fog Film

- Nosepiece and thin anti-fog film form smoothly across the nose and cheekbones while gently clinging to the skin to form a seal, preventing humid air from fogging up eyewear.

③ Smooth Naturally Breathable Inner Layer

- Provides a light and cool feel to the mask, reducing chances of irritation throughout extended wear.
- Inner layer maintains its structure, allowing the mouth to feel open and airy.



Not made with
natural rubber latex



* * * * *

VIVIT Supply LLC

200 Continental Dr., Suite 401
Newark, DE 19713 USA

W vivitsupply.com

V004