



MEDICAL FACE MASK

TIES • ANTI-FOG TAPE • GREEN

ASTM LEVEL 1



13313TM

For use in environments with minimal fluid exposure, this ASTM Level 1 mask establishes a physical barrier between clinicians and contaminants from their surrounding environment. The pressure-sensitive adhesive strip allows for a tailored, secure seal to prevent the fogging of eyewear, providing a clear line of sight during surgical procedures.

BFE	PFE	DELTA P $\Delta P(\text{mmH}_2\text{O}/\text{cm}^2)$	FLUID RESISTANCE	FLAME SPREAD
> 95%	> 80%	< 5.0	80 mm Hg	Class 1
MASKS/BOX		BOXES/CASE		MASKS/CASE
50		6		300

- ① **Premium Tie Material and Design**
 - Tie material that is comfortable on the skin throughout extended wear.
 - Long ties create a flexible fit while tying over top of a bouffant or surgical cap.
- ② **High Quality Nosepiece with Anti-fog Tape**
 - Nosepiece combined with a skin-friendly adhesive strip creates the ultimate anti-fog face mask, ensuring a complete seal across the nose and cheekbones to prevent the fogging of eyewear.
- ③ **Smooth Naturally Breathable Inner Layer**
 - Provides a light and cool feel to the mask, reducing chances of irritation throughout extended wear.
 - Inner layer maintains its structure, allowing the mouth to feel open and airy.

