

MEDICAL FACE MASK

EARLOOPS • BLUE

ASTM LEVEL 2



2011EM

Combining comfort and breathability, this ASTM Level 2 mask offers standard fluid resistance and filtration efficiency to keep clinicians protected during procedures and patient care. This mask can also be worn when experiencing respiratory symptoms, such as coughing or sneezing, to help limit the spread of droplets containing pathogens.

BFE	PFE	DELTA P $\Delta P(\text{mmH}_2\text{O}/\text{cm}^2)$	FLUID RESISTANCE	FLAME SPREAD
> 98%	> 85%	< 6.0	120 mm Hg	Class 1
MASKS/BOX		BOXES/CASE		MASKS/CASE
50		10		500

① Ultra-soft Earloop

- Earloop material is soft and smooth for maximum comfort.
- Round design makes this mask one of the most comfortable options available.

② High Quality Nosepiece

- Aluminum nosepiece forms smoothly across the nose and cheekbones.
- Remains securely in place and reduces fogging of eyewear.

③ Responsive Inner Layer Material

- Padded feel for ultimate comfort.
- Flexible material that moves with you.



VIVIT Supply LLC

200 Continental Dr., Suite 401
Newark, DE 19713 USA

W vivitsupply.com

V002