



MEDICAL FACE MASK

TIES • DUCKBILL • BLUE

ASTM LEVEL 2



20131TM

This ASTM Level 2 mask is designed for superior comfort and breathability, making it ideal for surgical procedures and extended wear. Its unique duckbill shape creates a large breathing chamber, keeping the mask away from the nose and mouth for easier airflow. Made from lightweight materials, it provides maximum comfort without compromising protection.

BFE	PFE	DELTA P ΔP(mmH ₂ O/cm ²)	FLUID RESISTANCE	FLAME SPREAD
> 98%	> 85%	< 6.0	120 mm Hg	Class 1
MASKS/BOX		BOXES/CASE		MASKS/CASE
50		6		300

① Premium Tie Material and Design

- Tie material that is comfortable on the skin throughout extended wear.
- Long ties create a flexible fit while tying over top of a bouffant or surgical cap.

② High Quality Nosepiece

- Malleable nosepiece forms smoothly across the nose and cheekbones.
- Remains securely in place and reduces fogging of eyewear.

③ Optimized Duckbill Chamber

- Large breathing chamber to keep the mask away from the nose and mouth during use.
- Provides maximum comfort, and breathability ideal for surgical procedures and extended wear.

