




MEDICAL FACE MASK HORIZONTAL TIES • BLUE ASTM LEVEL 2



20141TM

Designed for comfort and breathability, this Level 2 mask provides reliable fluid resistance and filtration efficiency to keep clinicians protected during procedures and patient care. Its wide fit enhances security with full coverage and a comfortable seal, while the horizontal tie design ensures a flexible fit with soft, skin-friendly ties for extended wear.

BFE	PFE	DELTA P $\Delta P(\text{mmH}_2\text{O}/\text{cm}^2)$	FLUID RESISTANCE	FLAME SPREAD
> 98%	> 85%	< 6.0	120 mm Hg	Class 1
MASKS/BOX		BOXES/CASE		MASKS/CASE
50		6		300

① Premium Tie Material and Design

- Tie material that is comfortable on the skin throughout extended wear.
- Long ties create a flexible fit while tying over top of a bouffant or surgical cap.

② High Quality Nosepiece

- Malleable nosepiece forms smoothly across the nose and cheekbones.
- Remains securely in place and reduces fogging of eyewear.

③ Wide Fit

- Extra coverage for wider faces, facial hair, or clinicians who prefer a larger mask.



Not made with
natural rubber latex



VIVIT Supply LLC

200 Continental Dr., Suite 401
Newark, DE 19713 USA

W vivitsupply.com

V001