



MEDICAL FACE MASK

EARLOOPS • OPAL

ASTM LEVEL 3



30114EM

Using only premium quality materials, this ASTM Level 3 mask offers the highest levels of fluid resistance and filtration efficiency. Rigorously tested against the ASTM F2100 Standards, clinicians can be confident they are protected during medical procedures and patient care in areas such as the ER and ICU.

BFE	PFE	DELTA P $\Delta P(\text{mmH}_2\text{O}/\text{cm}^2)$	FLUID RESISTANCE	FLAME SPREAD
> 98%	> 85%	< 6.0	160 mm Hg	Class 1
MASKS/BOX		BOXES/CASE		MASKS/CASE
50		10		500

① Ultra-soft Earloop

- Earloop material is soft and smooth for maximum comfort.
- Round design makes this mask one of the most comfortable options available.

② High Quality Nosepiece

- Aluminum nosepiece forms smoothly across the nose and cheekbones.
- Remains securely in place and reduces fogging of eyewear.

③ Smooth Naturally Breathable Inner Layer

- Provides a light and cool feel to the mask, reducing chances of irritation throughout extended wear.
- Inner layer maintains its structure, allowing the mouth to feel open and airy.



Not made with
natural rubber latex



* * * * *

VIVIT Supply LLC

200 Continental Dr., Suite 401
Newark, DE 19713 USA

W vivitsupply.com

V004