



# MEDICAL FACE MASK TIES • OPAL ASTM LEVEL 3



**30314TM**

Rigorously tested against the ASTM F2100 Standards, this ASTM Level 3 mask offers maximum protection for both fluid resistance and filtration efficiency. Clinicians can be confident when using this mask during surgical procedures including trauma, ER, and orthopedic cases.

BFE	PFE	DELTA P $\Delta P(\text{mmH}_2\text{O}/\text{cm}^2)$	FLUID RESISTANCE	FLAME SPREAD
> 98%	> 85%	< 6.0	160 mm Hg	Class 1
MASKS/BOX		BOXES/CASE		MASKS/CASE
50		6		300

## ① Premium Tie Material and Design

- Tie material that is comfortable on the skin throughout extended wear.
- Long ties create a flexible fit while tying over top of a bouffant or surgical cap.

## ② High Quality Nosepiece

- Aluminum nosepiece forms smoothly across the nose and cheekbones.
- Remains securely in place and reduces fogging of eyewear.

## ③ Smooth Naturally Breathable Inner Layer

- Provides a light and cool feel to the mask, reducing chances of irritation throughout extended wear.
- Inner layer maintains its structure, allowing the mouth to feel open and airy.



Not made with  
natural rubber latex



\* \* \* \* \*

**VIVIT Supply LLC**

200 Continental Dr., Suite 401  
Newark, DE 19713 USA

**W** [vivitsupply.com](https://vivitsupply.com)

V004