



MEDICAL FACE MASK

EARLOOPS • ANTI-GLARE / ANTI-FOG VISOR

ANTI-FOG FOAM • OPAL • ASTM LEVEL 3



311141EMV

Tested against the ASTM F2100 Standards, this ASTM Level 3 mask offers high filtration efficiency and is designed to withstand intense spray impact in areas such as trauma, ICU, and the ER. Its attached wrap-around visor with built in anti-fog properties and an anti-fog foam strip provides clinicians with optical clarity and enhanced protective features.

BFE	PFE	DELTA P $\Delta P(\text{mmH}_2\text{O}/\text{cm}^2)$	FLUID RESISTANCE	FLAME SPREAD
> 98%	> 85%	< 6.0	160 mm Hg	Class 1
MASKS/BOX		BOXES/CASE		MASKS/CASE
25		4		100

① Ultra Soft Earloop

- Earloop material is soft and smooth for maximum comfort.
- Round design makes this mask one of the most comfortable options available.

② High Quality Nosepiece with Ultra Soft Foam

- Aluminum nosepiece and soft foam strip work together to form smoothly across the nose and cheekbones.
- Remains securely in place to create a custom and secure seal providing enhanced resistance to fogging of eyewear and a cushion-like feel to the mask.

③ Anti-Fog Visor with Wrap-Around Design

- Visor design provides clinicians with a clear line of sight paired while protecting from direct and peripheral protection from spray of fluids.



Not made with
natural rubber latex



* * * * *

VIVIT Supply LLC

200 Continental Dr., Suite 401
Newark, DE 19713 USA

W vivitsupply.com

V002